


☐

I'm not robot


reCAPTCHA

I'm not robot!



Cobas c 311 operator manual

Cobas c311 operator manual pdf. Cobas c311 operator manual. Cobas c 311 analyzer. operator's manual.



roche diagnostics international ltd.

Table of contents 93 4 In this chapter 4 About positions of reagents



. 95 About statuses of reagents .



Analizador cobas c 311

Guía de seguridad – Versión 1.1

Versión del software 01-10



cobas

Light reader

. 98 Loading reagent packs 101 Unloading reagent packs

104 Changing ISE reagents

. 106 Replacing ISE reagents .

106 Performing a reagent prime .

. 108 Performing an ISE calibration

. 109 Changing sample probe detergents

112 Changing cell detergents

. 113 Performing a cell detergent prime 115 Defining reagent level alarms

. 116 Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 94 Table of contents Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Reagents 95 About positions of reagents You get an overview of all reagents used on the cobas c 311 analyzer and where to find them. Reagents in reagent packs A Gate for reagent packs B Reagent pack A B Reagents in reagent packs are used for photometric measurement, for special washing, and for maintenance.

They are located in the temperature-controlled reagent disk. You must load/unload reagent packs via the gate. Reagents in reagent packs are: • All photometric test reagents • Diluents (NaCl, A1CD2) • Detergents (NAOHD, SMS) • EcoTergent (mandatory reagent) Ecotergent is indicated as BathDet (water bath detergent) in the instrument software ISE reagents A B C ISE reagents are used for ISE measurement.

They are located in bottles between the reagent disk and the reaction disk. Reagents used for ISE measurement are: • ISE Internal Standard Gen.2 (ISE IS) • ISE Reference Electrolyte (ISE REF) • ISE Diluent Gen.2 (ISE DIL) A ISE DIL B ISE IS C ISE REF Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 96 About positions of reagents Evaporation of internal standard Downward drift in controls and patient sample results can occur due to increased evaporation of the internal standard. Evaporation of the internal standard can occur due to an unfavorable combination of: • Low humidity • High altitude • Air movement in the laboratory depending on the use of: • Air conditioning - Fans • Filling level of the internal standard In case of internal standard evaporation the use of a bottle insert for the internal standard is recommended.

For ordering the bottle inserts and information on the correct usage of the bottle inserts, please contact your Roche service representative. Calibrators A Sample container B CFAS A C ISE Standard (Low/High) B C Calibrators are prepared in sample containers and then placed on the sample disk. Controls used for photometric and ISE measurement are, for example: • PreciControl ClinChem 1 (PCCC 1) • PreciControl ClinChem 2 (PCCC 2) Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Reagents 97 ISE washing reagents A A Sample container B Activator C SysClean Sample probe detergents B C The ISE washing reagents are placed on defined positions in the sample disk. The Reagent used for cleaning the ISE flow path is: • ISE cleaning solution (SysClean) SysClean is placed on position W1. The Reagent used for conditioning the ISE electrodes is: • ISE Activator (Activator) Activator is placed on position W2. A B Sample probe detergents are located in bottles next to the sample disk. Detergents used for washing sample probes are: • Sample Cleaner 1 (SmpCln 1) (Multiclean) • Sample Cleaner 2 (SmpCln 2) (SMS) A SmpCln 1 B SmpCln 2 Cell detergents A A CellCln 1 B B CellCln 2 Cell detergents are located in bottles behind the front door of the analyzer. Detergents used for washing cells are: • Cell Wash Solution I / NAOHD (CellCln 1) • Cell Wash Solution II / Acid Wash (CellCln 2) u Related topics • Loading reagent packs (101) • Unloading reagent packs (104) • Replacing ISE reagents (106) • Changing sample probe detergents (112) • Changing cell detergents (113) Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 98 About statuses of reagents About statuses of reagents To ensure measurement, the required reagents must be on board and their fill level must be high enough.

Additionally, the reagent packs and the ISE reagents must have a valid calibration and QC. The Reagent Preparing button and the Reagent Overview button indicate if a reagent is empty, missing, or below a set volume level by changing its color. u Defining reagent level alarms (116) From the System Overview menu, you can choose the Reagent Overview button. q To view the Reagent Overview button in the System Overview menu, you must choose AU or ISE in the Overview group box. The Reagent Overview dialog box displays the statuses of all reagents on board. You can check the remaining amounts, the calibration statuses, and the QC statuses. Reagent packs The Reagent group box displays the reagent disk with the 42 reagent packs and their statuses. In the center of the reagent disk, you can view detailed information on the selected reagent pack. If you want to view more detailed information, you can choose the Detail button. A B C A Reagent pack B Detailed information C Color legend A C A Reagent pack status C QC status B Calibration status Each reagent pack is divided into three sections displaying the following statuses: • Reagent pack status • Calibration status • QC status Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Reagents 99 Section Status Cas Reagent status Cal QC Calibration status Quality control status Color The color of each section displays the current status. Description Active Reagent currently in use. Warning Reagent remaining tests below yellow alarm level. Zero tests Reagent empty. Reagent available in another reagent pack. Caution Standby Reagent empty.

Reagent not currently in use. Free Pos. Reagent position empty. Expired Valid Expiration date exceeded. Calibration valid. Requested Calibration requested. Failed Calibration failed. Valid Quality control valid. Requested Quality control requested. Violated Quality control failed. N.A. Quality control not applicable (no values available). y ISE reagents The ISE group box displays the fill level and the statuses of the ISE reagents ISE IS, ISE DIL, and ISE REF. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 100 About statuses of reagents A B C D Each ISE reagent is divided into four sections displaying the following statuses: • Fill level (remaining volume) as graphic display • Fill level (remaining volume) as amount (mL) • Calibration status • QC status A Fill level display C Calibration status B Fill level amount (mL) D QC status Section Status Cal Calibration status Color The colors of the calibration and QC sections display the actual statuses. Description Valid Calibration valid. QC Quality control status Requested Calibration requested. Failed Calibration failed. Valid Quality control valid. Requested Quality control requested. Violated N.A. Quality control failed. Quality control not applicable (no values available). y Detergents The Inventory group box displays the fill level of the detergents. A Fill level display A B B Fill level amount (mL) Each detergent is divided into two sections displaying the following statuses: • Fill level (remaining volume) as graphic display • Fill level (remaining volume) as amount (mL) Roche Diagnostics cobas c 311 analyzer

Software version 01-10 · Operator's Manual · 3.1 Reagents 101 Loading reagent packs If the remaining volume of a reagent pack is low, or the expiration date has expired, you must replace it. To ensure measurement, you must scan and insert a new reagent pack, or a reagent pack has expired. This analyzer scans the barcode label on the reagent pack. f The analyzer checks if the corresponding test application is available.

Reagents 102 Loading reagent packs If you insert a reagent pack the wrong way around, the reagent is pipetted in the wrong order, leading to incorrect results or system damage. r Always load the reagent packs in the direction shown on the label near the loading port. l WARNING Incorrect results due to insufficient reagent volume When using a cobas c pack MULTI, it must be filled with the exact filling volume specified in the package insert. If this is not the case, measurement accuracy may deteriorate. Additionally, an alarm (reagent short) may occur if the analyzer detects that the reagent volume is low. r Ensure that the cobas c pack MULTI contains the correct filling volume. Refer to the package insert. The analyzer monitors the remaining volumes by counting down. When loading new reagent packs, the analyzer starts counting at the initial number of available tests. When loading reagent packs that have already been used, the analyzer continues counting at the last monitored number of available tests. If a reagent is empty, missing or below a set volume level, the Reagent Preparing button and the Reagent Overview button change their color to red, yellow, or purple. u Defining reagent level alarms (116) You can print a Reagent Load/Unload List report that shows the reagent packs to be loaded. u Preparing the reagents (79) j The analyzer is in Standby status. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 102 Loading reagent packs 3 2 1 r To load reagent packs 1 Choose Reagent > Setting . 2 Press Ctrl and choose the reagents you want to load from the list. 3 Choose the Loading button. 4 5 4 Choose the Execute button. 5 Wait 1-2 minutes for the gate cover of the reagent disk to unlock. f You can see and hear the locking bar retract. 6 6 Open the gate cover. 8 7 7 With the reagent barcode label facing the barcode reader, put the reagent pack on the guide rail in front of the barcode reader. 8 To scan the reagent pack, slide the reagent pack along the guide rail to the rear stop. f The analyzer scans the barcode label on the reagent pack. f The analyzer checks if the corresponding test application is available.

Reagents 103 Reagents 103 6 6 Open the gate cover. The following data are encoded in this barcode: o System ID o Lot number o Cassette number o Expiration date o Bottle configuration information (reagent packs) 9 Wait for the reagent pack to be accepted. 10 With the barcode on the right side, load the reagent pack into the slot. 11 11 Close the gate cover. f The gate cover is locked automatically. 12 Wait up to 1 minute for the reagent pack to be registered.

If the reagent pack has not been on the analyzer before, the following actions are taken: • The analyzer registers the initial capacity. • The bottles of the reagent pack are pierced. f When the cassette is registered, you can see the cassette information in the Cassette Loading dialog box.

13 To finish loading, choose the End button. u Related topics • Defining reagent level alarms (116) Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 104 Unloading reagent packs Unloading reagent packs If you need an empty position on the reagent disk for loading a new reagent pack, or if the reagent is empty or below a set volume level, you must unload a reagent pack. In Reagent > Setting , in the Remaining Tests column, you can check for empty reagent packs. You can print a Reagent Load/Unload List report that shows the reagent packs to be unloaded. u Preparing the reagents (79) j The analyzer is in Standby status. 1 r To unload reagent packs 1 Choose Reagent > Setting .

2 Press Ctrl and choose the reagent packs you want to unload from the list. 3 Choose the Unloading button. 2 3 4 5 4 Choose the Execute button. 5 Wait up to 1 minute for the gate cover of the reagent disk to unlock. f You can see and hear the locking bar retract. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 7 Reagents 105 6 6 Open the gate cover. The following data are encoded in this barcode: o System ID o Lot number o Cassette number o Expiration date o Bottle configuration information (reagent packs) 9 Wait for the reagent pack to be accepted. 10 With the barcode on the right side, load the reagent pack into the slot. 11 11 Close the gate cover. f The gate cover is locked automatically. 12 Wait up to 1 minute for the reagent pack to be registered.

10 To finish unloading, choose the End button. 11 Store the reagent packs in a refrigerated place or dispose them according to local regulations. 1 You can reload a used reagent pack onto the same analyzer at a later time.

However, after deleting a test, you cannot reload used reagent packs for this test. Note: You must delete the test and download it again to, for example, change the unit of measure.

Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 106 Changing ISE reagents Changing ISE reagents ISE REF ISE DIL If the remaining volume of an ISE reagent is low or the expiration date has expired, you must replace it. When replacing the ISE reagent ISE IS or ISE REF, you must perform a reagent prime. Before starting operation, you must calibrate the ISE measuring system for all replaced ISE reagents. 1 Replace bottle 2 Prime 3 Calibrate In this section Replacing ISE reagents (106) Performing a reagent prime (108) Performing an ISE calibration (109) Replacing ISE reagents If the remaining volume of the ISE reagent ISE IS, ISE REF, or ISE DIL is low, or the expiration date has expired, you must replace it with a new and full ISE reagent bottle to ensure measurement. After replacing the ISE reagent ISE IS, ISE REF, or ISE DIL, you must reset the reagent volume. If you replace ISE IS or ISE DIL, the liquid level is detected automatically before the first measurement. If you replace ISE REF, the reagent volume is monitored by countdown. l WARNING Incorrect results due to ISE measurement failure r Perform ISE prime after resetting the ISE reagent volume. In the Reagent Overview dialog box, you can check the remaining volume of ISE reagents. You can print a Reagent Loading List report that shows the ISE reagents to be replaced. u About statuses of reagents (98) Preparing the reagents (79) j m The analyzer is in Maintenance mode. u Switching the analyzer to Maintenance mode (236) Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Reagents 107 1 r To replace ISE reagents 1 If necessary, insert the bottle insert into the new internal standard bottle. A B C 2 2 CAUTION! Incorrect result due to incorrect placement of reagents. Make sure to place the reagent bottles according to their correct positions. 2 Replace the corresponding bottle of ISE reagent with a new and full bottle.

A ISE DIL B ISE IS C ISE REF A B A Aspiration filter B ISE REF 3 3 If you replace the ISE REF bottle, clean the aspiration filter A in the tubing: • Take out the tubing from the bottle. • Unscrew the filter from the tube end. • Wash the filter with deionized water. • Then rinse the filter with deionized water. • Shake off the water. • Place the tubing in the new bottle, so that the end of the tube touches the bottom of the bottle. 4 Terminate Maintenance mode. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 108 Changing ISE reagents 6 5 5 Choose Reagent > Status . 6 On the Status tab, choose the changed ISE reagents. 7 Choose the ISE Volume Reset button. 7 8 To reset the fill volume to the level of a full bottle, choose the YES button. f The volume of ISE REF is displayed immediately. The volumes of ISE DIL and ISE IS are displayed with the next measurement. u Related topics • Cleaning the ISE REF aspiration filter (284) • Defining reagent level alarms (116) • Terminating the Maintenance mode (237) Performing a reagent prime When the tubing is disconnected from the ISE REF bottle, air may enter the electrodes. If the ISE IS bottle or the ISE REF bottle is replaced, the flow path must be filled with new liquid. In both cases, you must perform a reagent prime. n IS+REF: approximately 10 minutes • REF: approximately 2 minutes d m IS+REF, ISE IS 9230 u L, ISE REF 1690 u L, m REF, ISE REF 3120 u L j The analyzer is in Standby status. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Reagents 109 2 1 r To perform an ISE reagent prime 1 Choose Utility > Maintenance .

2 Select the (7) Reagent Prime option. 3 Choose the Select button. 3 4 5 7 4 Choose the Parameter button. 5 Select one of the following options: • To prime both the ISE measuring tubing and the ISE reference tubing, select the IS+REF option. • To prime the ISE measuring tubing, select the REF option. 6 Choose the OK button. 7 Choose the Execute button. f The reagent prime is completed when the analyzer returns to Standby status. Performing an ISE calibration If you replace ISE IS, ISE REF, or ISE DIL, you must calibrate the ISE.

According to your specific working requirements, you can define two different ISE calibration curves (e. g., type A for urine and type B for plasma). The ISE calibrators can be used with or without barcode. During ISE calibration, the ISE internal standard (ISE IS) is measured additionally. d m Calibrators ISE Low (S1), ISE High (S2), and ISE High (S3) ISE IS and ISE REF have been primed. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 110 Changing ISE reagents 2 3 1 r To perform an ISE calibration 1 Choose Calibration > Status . 2 From the Test column, choose ISE-A or ISE-B to select the calibration curve. 3 In the Method group box, choose the Full button. f The selected calibration method Full is displayed on the Calib. Method tab and highlighted green. f The Save button turns yellow. 4 4 Choose the Save button. f All tests with a highlighted green entry in the Calib. Method tab are requested for calibration.

5 Clean the filter in the tubing. • Remove the tubing from the cell detergent bottle. • Unscrew the filter from the tube end. • Wash the filter with deionized water. • Then rinse the filter with deionized water. • Screw the filter back on to the tube end. 3 3 Take a new and full cell detergent bottle. 4 Insert the tubing into the cell detergent bottle so that it touches the bottom of the bottle. 4 Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 114 Changing cell detergents 5 5 Position the new cell detergent bottle. 8 7 9 6 6 Choose Reagent > Status . 7 Choose the new cell detergent bottle. 8 Choose the Reagent Volume Reset button. f The Confirmation dialog box is displayed. 9 Choose the OK button. 10 Perform a cell detergent prime. u Related topics • Cleaning the detergent aspiration filters (286) • Performing a cell detergent prime (115) Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Reagents 115 Performing a cell detergent prime When the tubing is disconnected from the cell detergent bottles, air may enter the tubing. Therefore, you must perform a detergent prime. j The analyzer is in Standby status.

1 r To perform a cell detergent prime 1 Choose Utility > Maintenance . 2 Select the (8) Cell Detergent Prime option. 3 Choose the Select button. 2 3 5 4 4 Select the detergent to be primed. • Detergent 1 = CellCln 1 • Detergent 2 = CellCln 2 5 Choose the Execute button. f The detergent prime is completed when the analyzer returns to Standby status. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 116 Defining reagent level alarms Defining reagent level alarms To ensure that there is sufficient reagent on board for measuring, you can define reagent levels at which an alarm is issued. Color is alarm Meaning A reagent is empty. The remaining tests of a reagent are below the defined reagent caution level. The remaining tests of a reagent are below the reagent level defined for a purple alarm. y Reagent level alarms If an alarm is issued, in the System Overview menu, the Reagent Preparing button and the Reagent Overview button change the color. q To view the Reagent Overview button in the System Overview menu, choose the Reagent Overview button. q To view the Reagent Overview button in the System Overview menu, choose the Reagent Overview button.

Yellow alarm Purple alarm The red alarm is issued automatically. For yellow and purple alarms, you can define the alarm level for each reagent separately. • For test reagents, the alarm level refers to the remaining number of tests. • For ISE reagents, wash solutions, and diluents the alarm level refers to the remaining volume in mL. A yellow alarm is only issued if the reagent level alarm is defined for this test reagent, wash solution, or diluent.

A purple alarm is only issued under the following conditions: • The reagent level alarm has been defined for this test. • On the System Overview menu, the Preventive Action check box has been selected. It is recommended, that the purple alarm corresponds to the volume of reagents used in one day.

Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 3 Reagents 117 1 r To define a yellow reagent level alarm 1 Choose Utility > System 2 Choose the Page button until Page 2/5 is displayed on the button. 3 Choose the Reagent Level Check button. 2 4 6 8 4 In the Yellow Alarm group box, you can filter by reagent type: • The Assay option displays only test reagents. • The Without Assay option displays only wash solutions and diluents. 5 Choose the reagent for which you want to define a reagent level alarm. f The current alarm values are displayed in the Count column. 6 To activate the reagent level alarm for this test, select the Activate Alarm option. f The Yellow Check Update button turns yellow.

To confirm, choose the OK button. The new value in the Activate Alarm group box. • For the reagent packs, the level alarm lies between 10-500 tests. • For ISE reagents, wash solutions, and diluents, the level alarm lies between 5-9999 mL. 8 To save the changes, choose the Yellow Check Update button. 9 To confirm, choose the OK button.

Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 118 Defining reagent level alarms 3 1 r To define a purple reagent level alarm 1 Choose Utility > System 2 Choose the Page button until Page 2/5 is displayed on the button. 3 Choose the Reagent Level Check button. 2 4 6 8 4 In the Purple Alarm group box, you can filter by reagent type: • The Assay option displays only test reagents. • The Without Assay option displays only wash solutions and diluents.

5 Choose the reagent for which you want to define a reagent level alarm. f The current alarm values are displayed in the Count column. 6 To activate the reagent level alarm for this test, select the Activate Alarm option. f The Purple Check Update button turns yellow. 7 To change the alarm value of this test, enter the new value in the Activate Alarm group box. • For the reagent packs, the level alarm lies between 10-500 tests. • For ISE reagents, wash solutions, and diluents, the level alarm lies between 5-9999 mL. 8 To save the changes, choose the Purple Check Update button. 9 To confirm, choose the OK button. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 2 Publication information Publication version 1.0 1.1 2.0 3.0 3.0.1 3.0.2 3.1 y Revision history Software version Revision date 01-01 01-04 01-06 01-09 2007-11 2009-11 2010-12 2014-10 01-09 2016-07 01-09 01-10 2016-10 2016-12 Change description First version Addendum 3.0 implemented in Online Help General update Documentation update Change weight (kg) R-134a 0.102 CO 2 equivalent (tonne) 0.14 y Fluorinated greenhouse gas detail Global warming potential 1430 5 Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 6 Contact addresses Inside the European Union and EFTA member states Manufacturer of cobas c 311 instrument Authorized representative Hitachi High Technologies Corporation 1-24-14 Nishi-Shimbashi Minato-ku Tokyo 105-8717 Japan Roche Diagnostics GmbH Sandhofer Strasse 116 68305 Mannheim Germany Outside the European Union and EFTA member states Manufacturer of Hitachi High Technologies Corporation Roche Diagnostics GmbH Sandhofer Strasse 116 68305 Mannheim Germany Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Table of contents 7 Table of contents Publication information Contact addresses Table of contents Intended use Symbols and abbreviations What is new in publication version 3.1? Operation 1 Overview of the instrument About the c 311 analyzer Overview of the hardware Overview of the power switches 2 Overview of operation Overview of the user interface Using the online help Overview of cobas link 3 Before operation Overview of the main workflow Quick start guide Starting from Sleep mode Starting from power-off Troubleshooting system alarms Preparing the analyzer for operation Measuring patient samples 4 Reagents About positions of reagents About statuses of reagents Loading reagent packs Unloading reagent packs Changing ISE reagents Changing sample probe detergents Changing cell detergents Performing a cell detergent prime Defining reagent level alarms 5 Calibration About calibration About calibration methods About recommended calibration About manual calibration Performing recommended calibration Performing manual calibration Checking calibration results Assigning non-barcoded calibrators to sample disk positions 121 122 124 126 127 129 132 133 95 98 101 104 106 112 113 115 116 75 76 88 57 58 63 70 9 11 7 9 17 9 2 6 29 42 47 17 18 25 6 Quality control About quality control (QC) Performing QC for individual tests Performing QC after calibration Performing QC for individual standby 137 141 144 reagents 146 Checking QC results 149 150 About the handling of QC results Workflows for handling and deleting QC results Assigning non-barcoded controls to disk 152 positions 159 7 Orders and results Workflow for sample processing Ordering tests Loading samples Processing samples Working with sample results Processing returns Processing open test requests Unloading samples About barcodes Preparing a DVD disk 8 After operation Preparing the analyzer for switching off Switching off the analyzer 221 223 163 165 170 175 183 194 199 210 212 217 Maintenance 9 Overview of maintenance Maintenance definitions About maintenance pipes and maintenance pipe functions Checking maintenance status Requesting the Maintenance Report List of maintenance items 235 236 List of maintenance checks List of tools, materials, and solutions List of spare parts and replacement intervals 255 239 247 248 249 251 253 254 10 Periodic maintenance List of maintenance tasks Daily maintenance Weekly maintenance Monthly maintenance Every three months maintenance Every six months maintenance 259 261 270 278 295 301 309 Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 8 Table of contents 11 Maintenance as required Checking and draining the vacuum tank Replacing the sample probe Eliminating clogging of the sample probe Replacing the reagent probe 329 331 338 340 Eliminating clogging of the reagent probe Replacing nozzle tips on cell rinse nozzles 346 348 Eliminating clogging of the cell rinse nozzles 351 Appendix 12 Specifications System specifications Analytical specifications List of sample containers List of barcode specifications List of accessories and consumables 357 360 365 366 368 13 Glossary Index 371 Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Intended use The cobas c 311 analyzer is a fully automated, discrete clinical chemistry analyzer intended for the in-vitro quantitative and qualitative determination of analytes in body fluids. The cobas c 311 analyzer is intended to be used by trained laboratory technicians. The operational environment for the cobas c 311 analyzer are clinical laboratories, hospital laboratories and commercial hospitals, as well as private laboratories. Symbols and abbreviations Product names Symbols used in the publication Except where the context clearly indicated otherwise, the following product names and descriptions are included in this publication: Roche Diagnostics cobas c 311 analyzer cobas c pack y Product names Descriptor analyzer reagent pack n d j u r u f c p w o u Symbol t y z Explanation List item Related topics containing further information Tip Extra information on correct use or useful hints. Start of a task Extra information within a task Result of a user action within a task Frequency of a task Duration of a task Materials that are required for a task Prerequisites of a task Topic. Used in cross-references to topics. Task. Used in cross-references to tasks. Figure. Used in figure titles and cross-references to figures. Table. Used in table titles and cross-references to tables. Equation. Used in cross-references to equations. y Symbols used in the publication Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 9 10 Symbols used on products Symbol Explanation Global Trade Item Number Quantity contained in the package Quantity contained in the package Orientation of the package during transportation Abbreviations Complies with the directive 2011/65/EU on RoHS CSV EC EFTA EN IEC IVD LIS n/a QC SBS SD WEEE y Symbols used on products The following abbreviations are used. Abbreviation ANSI C0BI Definition American National Standards Institute Compendium of Background Information comma separated values European Community European Free Trade Association European standard International Electrical Commission In vitro diagnostic Laboratory information system not applicable Quality control scan before sample stop Standard deviation Waste Electrical and Electronic Equipment y Abbreviations Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 What is new in publication version 3.1? Filter function for Reagent Level Check Improved Review By Exception data alarm setting Data alarms for serum index checks Improved handling of invalid reagent packs Chapter structure From software version 01-10 on, you can filter the reagents by type on the Reagent Level Check window. The filter function allows you to only display diluents and wash solutions. This function simplifies setting the reagent level check values for the yellow and purple alarm. In addition, separate reagent level check values are now available for each reagent and each solution. u Defining reagent level alarms (116) On Utility > System > Page 2/5 > Review By Exception menu, you can select the reagents to be displayed. In version 01-10, the software has been improved so that all data alarms are now available on the Review By Exception window. Seven new data alarms have been implemented in software version 01-10. These data alarms provide more detailed information about a possibly problematic serum condition: lipemic, hemolytic, or icteric samples, or a combination thereof. Former software versions attached only the generic >Index data alarm to results generated from a lipemic, hemolytic, or icteric sample. In some cases, you must delete an application and download it again, for example, to change the unit of measure. You can no longer use reagent packs for this application that were already in use on the system before deleting the application. u Unloading reagent packs (104) The chapter structure was slightly adapted to be consistent with other Roche publications. For example, the chapter Specifications was moved to Part Appendix . u Specifications (355) 11 Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 12 Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Operation 1 Overview of the instrument 152 2 Overview of operation 273 3 Before operation 55 4 Reagents 93 5 Calibration 119 6 Quality control 135 7 Orders and results 161 8 After operation 219 Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Overview of the instrument Table of contents 15 1 Read this chapter to get familiar with the hardware of the instrument. In this chapter 1 About the c 311 analyzer. 17 Overview of the hardware 18 Top view 18 Front view 18 20 Front view with open front door 21 Rear view 23 Overview of the power switches. 25 Main circuit breaker. 25 Power and maintenance 26 Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 16 Table of contents Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Overview of the instrument 17 About the c 311 analyzer The cobas c 311 analyzer is an automated system for qualitative and quantitative in vitro diagnostics. The ionselective electrode (ISE) and the photometric measurement unit provide a wide variety of tests.

The cobas c 311 analyzer consists of the following units: • Analytical unit • Control unit A Analytical unit A B Control unit B Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 18 Overview of the hardware Overview of the hardware The cobas c 311 analyzer is a standalone device. For a smooth operation, it is important to have access to all parts of the analyzer. In this section Top view (18) Front view (20) Front view with open front door (21) Rear view (23) Top view You can access most of the components of the analytical unit from the top when opening the top cover. A B Sampling area Reaction area A B C ISE area D Reagent area C D Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Overview of the instrument 19 Sampling area Reaction area ISE area Reagent area Samples are manually loaded to the sample disk. They are scanned or placed to a manually assigned position. That way, they are identifiable and easy to follow. From the sampling area, samples are piped to the reaction or ISE area. A photometric measuring system measures the absorbance of the reaction mixture in reaction cells on the reaction disk. An ion-selective electrode measuring system measures the electromotive force (EMF) in millivolts between the electrode in the diluted sample solution and the electrode in the reference solution. The electrodes are piped to the reaction cells. They are scanned and stored in a refrigerated compartment. From the reagent area, reagents are pipetted to the reaction area. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 20 Overview of the hardware Front view You can access most of the components of the control unit from the front.

A B C D E A Touchscreen monitor B Keyboard/mouse C Height adjusting lever Touchscreen monitor Keyboard/mouse Height adjusting lever Personal computer F D Personal computer E Front door F Cooling unit The touchscreen monitor displays information.

You can navigate through the software and start the analyzer functions. With the keyboard or the mouse, you can navigate through the software and select or enter information. If you lift the height adjusting lever, you can adjust the height of the touchscreen monitor. The personal computer monitors the functions and operational modes of the analyzer. The installed software has a graphical user interface (GUI) that allows you to control all analyzer functions. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Overview of the instrument 21 Printer (optional) Front door Cooling unit If a printer is connected, you can print the results. Some components of the analytical unit are located behind the front door. The cooling unit cools the reagent packs in the refrigeration compartment. Front view with open front door You have to open the front door to access some components of the analytical unit. F A B C G H D E A Reagent syringe B Sample syringe C ISE sipper syringe D Cell cleaner I Reagent syringe E Cell cleaner II F Alarm volume knob G Water tank H Drain tube of the vacuum tank The reagent syringe controls the pipetting of reagent to the reaction cells. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 22 Overview of the hardware Sample syringe ISE sipper syringe Cell cleaner I Cell cleaner II Alarm volume knob Water tank Drain tube of the vacuum tank The sample syringe controls the pipetting of sample to the reaction cells. The ISE sipper syringe controls the pipetting of diluted samples from the reaction cells and the pipetting of ISE reagents to the ISE unit. Cell cleaner I (CellCln I) contains alkaline detergent (NaOH-D) for the cell rinse process. Cell cleaner II (CellCln II) contains acid detergent (acid wash) for the cell rinse process. You can adjust the sound volume of the alarm buzzer with the alarm volume knob. Deionized water supplies the cell rinse unit, the rinse stations, and the incubator bath. The reaction cells are emptied with vacuum: • The reaction mixture is aspirated and emptied into the highly concentrated waste container. • The reaction cell rinse water is aspirated and drained through the diluted waste line. With the drain tube, you drain remaining liquid from the vacuum tank. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Overview of the instrument 23 Rear view Some components are located at the rear side of the analyzer. A B C D E A Cooling fan B Power supply outlets (US only, behind cover) C Ethernet connector D Mains connection E Main circuit breaker Cooling fan Power supply outlets (US only) Ethernet connector Mains connection Main circuit breaker F F Water supply tube G ISE drain port H Waste line I Waste container J Condensation water tray The cooling fan cools down the analytical unit. The power supply outlets energize the control unit in case the power cable is too short to reach an external power source.

They are located behind a cover. There are two power supply outlets for the control unit. One is for the PC, and the other is for the monitor. The Ethernet connection allows communication between the analytical and the control unit. The mains connection energizes the analytical unit. The main circuit breaker powers the analyzer on and off. You should keep it switching the cool the contents. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 24 Overview of the hardware Water supply tube ISE drain port Waste line Waste container Condensation water tray The water supply tube supplies deionized water for diluting samples and for washing. This ISE drain port discharges highly concentrated waste. The waste line collects diluted liquid waste, for example, reaction bath liquid and detergents. The waste container collects highly concentrated liquid waste from the ISE unit.

A liquid level sensor controls the level of the liquid waste and issues an alarm if the container gets full. The condensation water tray collects water condensing at the cooling unit. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Overview of the instrument 25 Overview of the power switches You can switch on and off the entire analyzer or parts of the analyzer. In this section Main circuit breaker (25) Power and maintenance switches (26) Main circuit breaker With the main circuit breaker, you switch on and off the entire analyzer, including the reagent cooling unit. A A Main circuit breaker Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 26 Overview of the power switches Power and maintenance switches With the power and maintenance switches, you switch on and off special areas of the analyzer. A B C A Power switch for the analytical unit B Maintenance switch Power switch for the analytical unit Maintenance switch Power switch for the control unit C Power switch for the control unit With the power switch for the analytical unit, you switch on and off the analytical unit except for the reagent cooling unit. With the maintenance switch, you switch on and off the low voltage of the pipettors, the sample disk, and the reaction disk. With the power switch for the control unit, you switch on and off the personal computer. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Overview of operation Table of contents 27 2 In this chapter 2 Overview of the user interface 37 About screen sections 29 About the System Overview menu 32 About the main menus 34 About status colors. 37 About buttons 38 List of analyzer statuses. 39 List of keyboard shortcuts. 40 Using the online help 42 Displaying the online help 42 Searching the online help 42 Searching via table of contents 43 Searching via index 43 Searching the full text 44 Full-text search options 44 Using favorites in the online help 45 Adding a favorite 45 Displaying a favorite 46 Overview of cobas link 47 About cobas link 48 About the cobas link information flow. 50 About cobas eLibrary 52 Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 28 Table of contents Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Overview of operation 29 Overview of the user interface To operate the analyzer, you must be familiar with the user interface. In this section About screen sections (29) About the System Overview menu (32) About the main menus (34) About status colors (37) About buttons (38) List of analyzer statuses (39) List of keyboard shortcuts (40) About screen sections A D B A B C Status line Menu area C Help line D Global buttons Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 30 Overview of the user interface Menu area - System Overview The user interface is divided into four sections: • Menu area • Status line • Help line • Global buttons The menu area displays the System Overview menu or one of the five main menus. To operate the user interface, use the touchscreen, mouse, or keyboard. The System Overview menu is the root menu of the user interface. You start daily work from this menu. u About the System Overview menu (32) Menu area - main menu In the five main menus you enter, update, and evaluate data. u About the main menus (34) Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Overview of operation 31 Status line The status line contains the following information: • Status of host, core unit, and analytical unit • Analyzer status • Logon name • Date and time When performing a maintenance task, the status line also displays the remaining time. Help line The color-coded buttons in the status line are only displayed with the main menus. The color of each button indicates the status of the respective module.

Indices status of the analytical unit. The LIS is displayed only when an LIS is connected. About status colors (37) List of analyzer statuses (39) The help line contains the Help button and a hint. To display the online help, choose the Help button. The online help is context-sensitive and contains the full Operator's Manual and the Software Description. Hints depend on the location of the cursor and suggest actions or types of information to enter.

Using the online help (42) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 32 Overview of the user interface Global buttons Global buttons provide quick access to key functions of the analyzer. Stop Shutdown S.Stop Alarm Media Eject Print Pause/Scan Start To stop operation in case of emergency. To log off, start Sleep mode, or shut down the analyzer. To stop sampling. The analyzer completes measuring already pipetted samples. To display the list of current alarms. If an alarm was issued, the Alarm button flashes yellow or red. To disconnect a USB flash drive. To print predefined lists and reports. To pause sampling, for example, to insert a STAT sample. To scan the sample containers on the sample disk. To start sampling. u Related topics • About the System Overview menu (32) • About the main menus (34) • About status colors (37) • List of analyzer statuses (39) • Using the online help (42) About the System Overview menu The System Overview menu is the root menu of the user interface. You start daily work from this menu. System Overview You use the System Overview menu to prepare the analyzer for daily operation. The menu also displays information about the modules of the analyzer and provides access to the five main menus. u About the main menus (34) Quick start guide (58) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Overview of operation 33 Work Flow Guide The Work Flow Guide group box supports daily routine tasks. From left to right, six buttons represent tasks required during daily operation. If a button turns yellow or red, you must act. To display the dialog box for the required action, choose the respective button. u About status colors (37) Overview group box The Overview group box contains color-coded status information on each module of the analyzer. To display status information, choose one of the following buttons: Host To display the status of the LIS (if connected), Core Control Unit AU To display the status of the control unit, To display the status of the analytical unit. ISE Alarm To display the status of the ISE unit. To display the Alarm dialog box. Reagent Overview To display the Reagent Overview dialog box. Main menu buttons To display a main menu, choose one of the following buttons: Workplace Reagent Calibration QC Utility To display the Workplace menu. To display the Reagent menu. To display the Calibration menu. To display the QC menu. To display the Utility menu.

Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 34 Overview of the user interface u Related topics • About screen sections (29) • About the main menus (34) • About status colors (37) • Quick start guide (58) About the main menus Main menu buttons In the five main menus, you enter, update, and evaluate data. To display one of the five main menus, choose one of the following buttons in the System Overview menu: • Workplace • Reagent • Calibration • QC • Utility Workplace menu The Workplace menu has the following submenus: Test Selection To order tests for samples. Data Review To review and manage test results. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Overview of operation 35 Reagent menu The Reagent menu has the following submenus: Setting Status To load and unload reagents. To check the number of available tests and to reset the volumes of sample cleaners, cell cleaners, and ISE reagent. Calibration menu The Calibration menu has the following submenus: Status Calibrator Install To review calibration status. To assign a non-barcoded calibrator to a position on the sample disk. To add a calibrator. To update calibrator data. QC menu The QC menu has the following submenus: Status Run Status Manual Currents To view the current status of the calibration. To view the calibration QC chart. To assign a non-barcoded control to a position on the sample disk. To add a control and to update QC data. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 36 Overview of the user interface Utility menu Return to System Overview The Utility menu has the following submenus: System Wash To configure the basic settings of the analyzer, for example, the use of sample disk positions. To page through the five pages of the submenu, choose the Page 1/5 button. Maintenance To select and monitor maintenance tasks. Application To download test applications from cobas link. Calc. Test To define formulas for calculated tests. Special Wash To program special wash instructions for the sample probe, reagent probe, and the reaction cells. Report Format To define the patient report form. Module Set To assign tests to the analyzer. If you have an operator level password only, you cannot edit all the settings. To return to the System Overview menu, choose one of the buttons in the status line. u Related topics • About screen sections (29) • About the System Overview menu (32) • About Maintenance mode (236) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Overview of operation 37 About status colors Status colors in System Overview To assure smooth operation of the analyzer, you must pay attention to the status colors. Status colors alert you to take timely action. Status colors are used in the Work Flow Guide group box and the Overview group box of the System Overview menu. Status colors in status line Status colors are also used in the status line above the main menus. Color coding The general meaning of the colors is as follows: Module in operation or maintenance mode. Module in standby mode. Action needed today (only for reagent preparing). Action needed soon. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 38 Overview of the user interface Action needed immediately. Examples The communication with the host is operational. The analytical unit is in standby mode. Action needed today (only for reagent preparing). Action needed soon. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 39 Overview of operation 39 Sort button In the column header, Sorts list by the column. Not all columns have this button. Examples Global button. Displays the Alarm dialog box. Red indicates that immediate action is needed. Saves changes in a dialog box or menu. Yellow indicates that you changed some settings and must now save those. Confirms a selection and closes the dialog box. Displays the Reagent Overview dialog box. Red indicates that immediate action is needed. Workflow button on the System Overview menu. Displays the Sample Tracking dialog box. Button in the status line. Displays the System Overview menu. Sort button in the header of the reagent, calibrator, and control status lists. Sorts these lists by test name. A white sort button marks the column that the list is sorted by. u Global buttons (32) List of analyzer statuses The analyzer status (mode) is displayed in the status line. The following list describes the analyzer statuses: Power Up The analyzer is loading programs, performing self-checks, etc. Initialize Standby After power-up, the analyzer initializes itself. The analyzer is ready for operation, but is not pipetting, measuring, or performing an automatic maintenance task. Preparation The analyzer is preparing for operation, for example, washing cuvettes. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 40 Overview of the user interface Operation Sample Reagent The analyzer is ready for operation. The analyzer is waiting for more samples. Sampling Pause Sampling Stop Stop The analyzer pauses sampling. You may load more samples on the sample disk. The analyzer has finished pipetting and is still measuring the samples. Standby (yellow background) The maintenance switch is on. The analyzer is in maintenance mode. You may open the top cover as all components on top are switched off. Shutdown The analyzer has performed an emergency stop due to hardware failure, or because any of the safety devices requested an emergency stop. The analyzer requires either a complete power-off or a reset (Utility > Maintenance) to resume operation. The analytical unit and the control unit are disconnected from the power supply. However, power for cooling the reagents is supplied. u Related topics • Status line (31) • About Maintenance mode (236) List of keyboard shortcuts F6 F7 F4 F5 F1 F2 F3 F8 F9 F10 F11 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 For quick access to frequently used menus and dialog boxes, use the following keyboard shortcuts: To display the online help. To display the Start Conditions dialog box. To display the Sample Stop dialog box. To display the Stop dialog box. To display the Workplace menu. To display the Reagent menu. To display the Calibration menu. To display the QC menu. To display the Utility menu. To display the Sample Tracking dialog box. To display the Alarm dialog box. Overview of operation 41 F12 To display the System Overview menu. Print Screen To display the Print dialog box. To take a screenshot, press Shift + Print Screen. Scroll Lock To display the Cancel Print dialog box. Pause/Break To display the Shutdown dialog box. Esc To exit a dialog box. u Related topics • Global buttons (32) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 42 Using the online help To find information quickly, use the online help. The online help comprises the full Operator's Manual and the Software Description . In this section Displaying the online help (42) Searching the online help (42) Using favorites in the online help (45) Displaying the online help 1 The online help is context-sensitive. It contains detailed information about the current menu or dialog box. r To display the online help 1 Navigate to a menu or dialog box you need information about, for example, choose Calibration > Status . 2 Choose the Help button or press F1 . 2 Searching the online help To find information in the online help, pick one of the following search methods. In this section Searching via table of contents (43) Searching via index (43) Searching the full text (44) Full-text search options (44) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Overview of operation 43 Searching via table of contents 1 To find a topic you know about, search via the table of contents. r To search via table of contents 1 In the online help, choose the Contents tab. 2 On the navigation pane, choose the desired chapter and topic. 2 Searching via index 3 2 1 To find a key word presumably included in the index, search via the index. r To search via index 1 In the online help, choose the Index tab. 2 In the text box, type a key word. f A list of matching index entries is displayed. 3 Choose the desired index entry. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 44 Using the online help Searching the full text 2 1 To find any word or phrase, search the full text. r To search the full text 1 In the online help, choose the Search tab. 2 In the text box, type the word or phrase you are looking for. 3 4 Full-text search options 3 Choose the List Topics button. f A list of matching topics is displayed. 4 Choose the desired topic. For an advanced full-text search, use the following logical search options: AND OR NEAR NOT To search for more than one word or phrase in a topic, use this operator. To search for any of two words or phrases in a topic, use this operator. To find words within ten words of each other, use this operator. To exclude topics with a particular word from the search, use this operator. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Overview of operation 45 Search previous results Match similar words To limit your search to the previous result list, select this check box. To find similar words and phrases, for example, plural forms, select this check box. Search titles only To search in titles only, select this check box. Using favorites in the online help To revisit topics easily, use favorites. In this section Adding a favorite (45) Displaying a favorite (46) Adding a favorite To create a list of favorites, add your favorite topics to the list. 1 r To add a favorite 1 In the online help, navigate to the topic you want to revisit later. 2 2 Choose the Favorites tab. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 46 Using the online help 3 3 If you want to rename the topic for the list of favorites, enter a new topic title in the text box. 1 To organize the list of favorites, use a prefix for topic titles (e.g., 1 Workplace topic X, 2 Reagent topic Y). 4 Choose the Add button. f A new favorite is added to the list of favorites. Displaying a favorite 1 r To display a favorite 1 In the online help, choose the Favorites tab. f The list of favorites is displayed. 2 Choose the desired favorite. 2 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Overview of operation 47 Overview of cobas link cobas link facilitates your work in several ways: • Automatic downloads of current applications, calibrators, and controls from the Roche remote service platform to your cobas instrument • Quick access to instructions for use, value sheets for calibrators and controls, etc. in cobas e-library • Electronic archiving of e-library package inserts and instrument settings • Automatic backup of instrument data to the cobas link data station • Upload of statistical data from the cobas instruments for performance monitoring and QC management • Screen sharing with a Roche Service representative for troubleshooting and support In this section About cobas link (48) About the cobas link information flow (50) About cobas e-library (52) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 48 Overview of cobas link About cobas link To access the Roche remote service platform, your cobas instruments must be connected to a cobas link. Remote service platform Remote service platform Internet cobas link c311 The Roche remote service platform provides information services for you and your cobas instruments via a secure Internet connection. The following applications are available: cobas e-library cobas e-support Automatic download of the most recent information on applications, calibrators, and controls. Screen sharing with a Roche Service representative for troubleshooting and support. cobas e-LabPerformance Upload of statistical data from the cobas instruments for performance monitoring and QC management. These applications are referred to as cobas e-services. cobas link Remote service platform Internet cobas link cobas link is a software package running on a dedicated computer in the laboratory, the cobas link data station. The data station and the software are together referred to as a cobas link. A cobas link connects one or more cobas instruments to the remote service platform. cobas links are set up and managed by the Roche Service representative. About cobas link (48) Overview of operation 49 Cobas link data station Cobas link data station The cobas link is connected to the control unit of the cobas c 311 instrument. c311 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Overview of operation 49 Remote service platform Offline mode cobas link in offline mode, the cobas link has no connection to the Internet. The cobas link cannot download information from the remote service platform or upload information. Software patches for cobas link are implemented via update DVDs. cobas e-library is updated via update CDs. The cobas link releases electronic barcodes to the instruments. The cobas link stores backups of the instrument data. c311 Remote service platform Online mode Internet cobas link In online mode, the cobas link has a stable connection to the Internet and connects the instruments in the laboratory to the remote service platform. The cobas link downloads information from the remote service platform and releases it to the instruments. The cobas link collects information from the instruments and uploads it to the remote service platform. u Related topics • About the cobas link information flow (50) • About cobas e-library (52) c311 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 50 Overview of cobas link About the cobas link information flow The information flow between the Roche remote service platform and the cobas link, and between the cobas link and the cobas instruments is bidirectional. Downloads from remote service platform Remote service platform Downloads from the remote service platform include the following: • cobas e-library packages: contain e-library package inserts and e-barcodes. • cobas link software updates Internet cobas link e-library package inserts e-barcodes The following documents come as e-library package inserts: • Instructions for use • Value sheets for calibrators and controls • Important notes, for example, about reassigned control values • Announcements from the local technical support You can read e-library package inserts in the cobas e-library application on the cobas link data station. u About cobas e-library (52) e-barcodes contain machine-readable information for cobas instruments: • Test-specific instrument settings • Lot-specific instrument settings • Values of calibrators and controls e-barcodes are stored on the cobas link data station. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Overview of operation 51 Parameter downloads from cobas link cobas link The instrument software indicates whether you must download new e-barcodes from the cobas link to the control unit of the instrument. c311 You must download e-barcodes if the Parameter Download button of the System Overview menu is red. u Downloading parameters from cobas link (86) Uploads to cobas link cobas link c311 A daily backup of the instrument data must be uploaded to the cobas link. • The backup is performed automatically if the cobas link Essential Information Upload maintenance item is included in the daily maintenance pipe. • cobas link stores the last five backups. • The backups are used to restore the control unit of the instrument, for example, in case of a hard disk crash. u List of maintenance items (251) u Related topics • About cobas link (48) • About cobas e-library (52) • Downloading parameters from cobas link (86) • List of maintenance items (251) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 52 Overview of cobas link About cobas e-library The cobas e-library application is an electronic archive running on the cobas link data station. You must check daily the list of new entries for the latest cobas e-library inserts. w cobas e-library user interface Access to cobas e-library Content of cobas e-library e-library packages You can access cobas e-library on the cobas link data station. For each application, you can download and control cobas e-library packages and the preceding version. Any older e-library packages are deleted from the cobas e-library. cobas e-library is updated daily via an automatic download from the Roche remote service platform. e-library packages contain e-library package inserts and e-barcodes. e-library package inserts are for your information. • e-barcodes contain machine-readable parameters and are downloaded to the cobas instruments. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Overview of operation 53 e-library package inserts The following documents come as e-library package inserts: • Instructions for use • Value sheets for calibrators and controls • Important notes, for example, about reassigned control values • Announcements from the local technical support List of new entries In the New Entries menu, the list of new entries is limited to 200 entries. Entries are not older than 30 days. You can filter the list, for example, for unread documents related to your cobas instrument. To find entries older than 30 days, you must search the cobas e-library archive. Other features of cobas e-library Other features of cobas e-library include the following: Comment Search Preferences You can comment on documents, for example, with your initials, to mark a document you have read. You can search for documents, for example, via the catalogue number on a reagent cassette. You can define preferences about documents that you want to see displayed in the New Entries menu. Online help To display the online help, choose the Help button on the bottom left of the cobas e-library screen. The online help is context-sensitive and contains the Operator's Manual on cobas e-library. u Related topics • About cobas link (48) • About the cobas link information flow (50) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 54 Overview of cobas link Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Before operation Table of contents 55 3 In this chapter 3 Overview of the main workflow 57 Quick start guide 58 Starting from Sleep mode 59 Starting from power-off 60 Performing maintenance and switching on the analyzer 61 Troubleshooting system alarms 62 Troubleshooting system alarms 63 Waking up the analyzer automatically 64 Waking up the analyzer manually 65 Monitoring the execution of the daily pipe 66 Checking the analyzer and performing maintenance 66 Starting from power-off 70 Performing maintenance and switching on the analyzer 71 Troubleshooting system alarms 72 Troubleshooting system alarms 73 Preparing the analyzer for operation 74 Preparing the analyzer for operation 75 Preparing the analyzer for operation 76 Performing daily maintenance 76 Clearing the sample database 78 Preparing the reagents 79 Performing recommended calibration and QC 81 Requesting recommended calibration and QC 82 Printing load lists for calibrators and controls 83 Measuring calibrators and controls 84 Downloading parameters from cobas link 86 Measuring patient samples 88 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 56 Table of contents Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Before operation 57 Overview of the main workflow The following chart summarizes the quick start guide. 22 19 Top (Standby mode) Measure patient samples Load/unload reagents Software FREDVFPDQZRUNARZ Top (Maintenance mode) Back Front/Slide 1 Enter login data 2 Check photometer report 3 Check syringes and tubing system 5 Check for water splashes 4 Switch to maintenance mode 21 Perform recommended calibration and QC 20 Download parameters from cobas link 6 Clean analyzer surface 7 Check reaction cell covers 13 Perform other due maintenance 8 Clean cell rinse nozzles 12 Empty liquid waste container 18 Clear sample database 17 Perform softwarecontrolled maintenance tasks 16 k&HFNZRUNARZ buttons 15 Check system alarms 9 Clean shield pipe 10 Clean probes and ISE sipper 11 Clean ISE drain port 14 Terminate maintenance mode u Related topics • Quick start guide (58) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 58 Quick start guide Quick start guide Step 1 Enter your login data. Announcements To prepare the analyzer for daily operation, you must perform the following steps. 1. Choose the Start Up button on the Reagent Load/Unload List. 3. Unload/load the reagents according to the report. If the Calibration and QC Select button is yellow, 1. Choose Calibration and QC Select > Recommended . 2. If a start-up group is defined, choose the Start Up button to request calibration for this group as well. 3. If the last shift started less than 24 hours ago, request the ISE calibration manually. 4. Choose the Routine QC button. 5. Choose the Stand By Bottle QC button. 6. Print the load lists for the recommended calibrators and controls. 7. Place the calibrators and controls on the sample disk. Use only fresh calibrators and controls. 8. Start the calibration and QC measuring process. 9. Check the calibration results. If the Parameter Download button is red, perform the following: 1. Choose the Parameter Download button. 2. Download the required parameters. 3. Repeat the calibration and QC measurements for the corresponding tests. 62 Quick start guide Step 22 Measure patient samples. y Quick start guide User action The following procedure assumes that test orders are sent from an LIS. 1. Deselect the Preventive Action check box. 2. Check that the green ACCESS SAMPLE DISK lamp is on. 3. Place the sample containers into the appropriate sample disk positions. 4. Start the measuring process. 5. To track sample processing, choose System Overview > Sample Tracking. 6. To review results, choose Workplace > Data Review. u Related topics • Starting from Sleep mode (63) • Starting from power-off (70) • Preparing the analyzer for operation (76) • Replacing reaction cells and cleaning incubator bath (278) • Replacing the photometer lamp (313) • Measuring patient samples (88) • About Maintenance mode (236) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Before operation 63 Starting from Sleep mode 1 Start analyzer from Sleep mode Automatic Automatic wake-up at set time and execution of daily pipe Manual 1 Wake-up analyzer manually 2 Monitor execution of daily pipe Typically the analyzer is in Sleep mode overnight. A wakeup time is set. In this case, the analyzer wakes up automatically and executes the daily pipe. After the daily pipe is completed, you must check the analyzer and perform the required manual maintenance tasks. However, you can also perform the manual maintenance tasks at the end of each work shift. If you wake up the analyzer manually before the set wake-up time, you additionally monitor the correct execution of the daily pipe. 2 Check analyzer and perform manual maintenance 3 Check analyzer and perform manual maintenance What is the daily pipe? The daily pipe is a maintenance pipe. A maintenance pipe executes a set of software-controlled maintenance items. The daily pipe includes all maintenance items required for instrument start. Typically the analyzer is set up to perform the daily pipe automatically at wake-up or power-up of the analyzer. If you do not use a daily pipe, you must manually execute the maintenance items included in the daily pipe. u Daily pipe (244) In this section Waking up the analyzer automatically (64) Waking up the analyzer manually (65) Monitoring the execution of the daily pipe (66) Checking the analyzer and performing maintenance (66) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 64 Starting from Sleep mode Waking up the analyzer automatically Typically the analyzer wakes up from Sleep mode automatically at the set wake-up time. To save time for routine operation, we recommend setting the wake-up time two hours before the work shift begins. The analyzer automatically performs all required maintenance items and is ready for operation when you arrive. j m The wake-up time is set. m The top cover is closed. 1 r To wake up the analyzer automatically 1 At the set wake-up time—while you are absent—the control unit PC performs a shutdown. f The shutdown deletes all temporary Windows files and the print buffer for reports. f The restart takes about 6 minutes. f After initialization, the daily pipe is executed. 2 2 When you arrive at the analyzer, enter your login data and choose the Logon button. f The System Overview menu is displayed. 3 3 Wait until Standby status is displayed. This status indicates that the daily pipe is complete. 4 Check the analyzer. u Related topics • Checking the analyzer (72) • About Shutdown status and Sleep mode (223) • List of recommended maintenance pipes (244) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Before operation 65 Waking up the analyzer manually If you want to start before the set time, you wake up the analyzer manually. c When you start before the set wake-up time j m The analyzer is in Sleep mode. M The top cover is closed. 1 r To wake up the analyzer manually 1 During Sleep mode, the Logon dialog box is displayed. Enter your login data and choose the Logon button. 2 2 Select the With Wake Up Pipe check box. f The control unit PC performs a shutdown. The shutdown deletes all temporary Windows files and the print buffer for reports. f The restart takes about 6 minutes. f After initialization, the daily pipe is executed. 3 3 When the Logon dialog box is displayed again, enter your login data and choose the Logon button. f The System Overview menu is displayed. 4 Monitor the execution of the daily pipe. u Related topics • About Shutdown status and Sleep mode (223) • List of recommended maintenance pipes (244) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 66 Starting from Sleep mode Monitoring the execution of the daily pipe If you start the instrument manually, you must monitor the execution of the daily pipe. r To monitor the execution of the daily pipe 1 Choose Utility > Maintenance . 1 2 2 Choose the Monitor button. 3 In the Maintenance Monitor dialog box, observe the execution of the maintenance items until (5) Air Purge is started. 4 Check the analyzer and perform manual maintenance. Checking the analyzer and performing maintenance After the daily pipe has been completed, you must check the functional integrity of the photometric components. You also check the analyzer surfaces and perform the daily maintenance tasks if this has not been done at the end of each work shift. d m Disposable cloth or paper towels m Laboratory disinfectant (no bleach) j The analyzer is in Standby status. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 2 Before operation 67 3 1 r To check the analyzer and perform manual maintenance 1 Choose Utility > Maintenance . 2 From the left list, choose the maintenance pipe programmed for wake-up or power-on (daily pipe). 3 To ensure that the daily pipe was executed on this day, check the displayed date. 4 4 Choose the Print button. 5 5 Choose the View button. 1 The (2) Photometer Check automatically generates a report. Therefore, you do not need to select the Photometer Check option. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 68 Starting from Sleep mode 6 7 6 In the Photometer Check report, check the current absorbance values. If the current results exceed 14 000 absorbance units at any wavelength, check the following points: • Make sure that the reaction cells, the incubator bath, and the photometric windows are free of contamination or bubbles. • Make sure that the reaction cells are at least half filled with water. • Then repeat the maintenance item (2) Photometer Check. If the results do not improve, replace the photometer lamp. 7 Open the front door. A A A Syringe connections B Syringe B 8 8 WARNING! Incorrect results due to inaccurate pipetting. Leaking syringe connection or air bubbles in a syringe may lead to incorrect results. 8 Check all syringes for leakages and air bubbles. • If a syringe connection is leaking, fasten it fingertight. • If there are air bubbles in the syringe, tap your finger lightly against the syringe while the liquid is flowing. • If you cannot remove the leakage or air bubbles, contact your Roche Service representative. 9 Make sure that no other tube connections or containers are leaking. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Before operation 69 10 10 Turn the maintenance switch to Maintenance mode. A A Maintenance switch A Rinse station of sample probe A B B Rinse station of reagent probe 11 13 11 Check the areas around the rinse stations A and B for water splashes. Water splashes may occur during the maintenance item (5) Air Purge. Splashes indicate that a probe is slightly clogged or worn. • If you detect splashes, ensure the analyzer is in Maintenance mode. • Wipe up any splashes with a paper towel moistened with disinfectant. • Eliminate clogging from the probe. 4 Open the front door. A B B 5 5 WARNING! Incorrect results due to inaccurate pipetting. Insert or remove a USB flash drive only in Standby status. • Use only one USB flash drive at a time. • Before removing a USB flash drive, choose the Media Eject button. • To prevent a virus from infecting the software, use the USB flash drive exclusively on the analyzer. • Wait until Standby status is displayed with white background. • Choose Utility > Maintenance and perform a reset. u Related topics • Switching the analyzer to Maintenance mode (236) • Eliminating clogging of the sample probe (338) • Eliminating clogging of the reagent probe (346) • Replacing the photometer lamp (313) • List of maintenance tasks (259) • Daily maintenance (261) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 70 Starting from power-off Starting from power-off Start analyzer from power off 1 Perform maintenance and switch on analyzer 2 Monitor execution of daily pipe 3 Check analyzer after completion of daily pipe If the analyzer is switched off at the beginning of your work shift, you must perform the manual maintenance tasks before you switch on the instrument. Identical to the start from Sleep mode, you monitor the execution of the daily pipe and check the functional integrity of the analyzer afterwards. In this section Performing maintenance and switching on the analyzer (70) Checking the analyzer (72) Performing maintenance and switching on the analyzer To ensure that the instrument surfaces are clean, you must perform the manual maintenance tasks. When switching on the analyzer, you must maintain the correct order of steps. If you perform manual maintenance as long as the analyzer is off, you do not need to switch to Maintenance mode. u Daily maintenance (261) j The analyzer is switched off. f To perform maintenance and switch on the analyzer 1 Perform all manual daily and other required maintenance tasks. 2 Make sure that the valve for the water supply is open. 2 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 3 3 Switch on the analyzer unit. 4 4 Switch on the analyzer unit The analyzer starts initialization. 5 5 Switch on the monitor and the printer. f After six minutes, the Logon dialog box is displayed. 6 6 Enter your login data and choose the Logon button. f After initialization, the analyzer starts the daily pipe by the Power Up Pipe function. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 72 Starting from power-off 7 7 Wait until Standby status is displayed. This status indicates that the daily pipe is complete. 8 Check the analyzer. u Related topics • Checking the analyzer (72) • List of maintenance tasks (259) • Daily maintenance (261) • Using the Power Up Pipe function (242) Checking the analyzer After the daily pipe has been completed, you must check the functional integrity of the photometric components with the Photometer Check report. You must also check the analyzer surfaces for water splashes. d m Disposable cloth or paper towels m Laboratory disinfectant (no bleach) r To check the analyzer 1 Choose the Print button. 1 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Before operation 73 2 2 Choose the View button. 1 The photometer check automatically generates a report. Therefore, you do not need to select the Photometer Check option. 3 3 In the Photometer Check report, check the current absorbance values. If the current results exceed 14 000 absorbance units at any wavelength, check the following points: f Make sure that the reaction cells, the incubator bath, and the photometric windows are free of contamination or bubbles. f Make sure that the reaction cells are at least half filled with water. f Then repeat the maintenance item (2) Photometer Check. If the results have not improved, replace the photometer lamp. • Open the front door. A B B 5 5 WARNING! Incorrect results due to inaccurate pipetting. Insert or remove a USB flash drive only in Standby status. • Use only one USB flash drive at a time. • Before removing a USB flash drive, choose the Media Eject button. • To prevent a virus from infecting the software, use the USB flash drive exclusively on the analyzer. • Do not store any other data on this USB flash drive. 1 r To back up the database before clearing the Sample Data Clear button. 2 Connect a USB flash drive to one of the USB ports or insert a DVD. 1 If a DVD and a USB flash drive are available, the data is stored on the USB flash drive. On a DVD, only 2 000 data sets can be stored. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Before operation 79 3 3 Select the Backup and Clear option and choose the OK button. Preparing the reagents 4 6 4 To save the data on a DVD or a USB flash drive, select the Mass Storage option. • Select a file name, for example, the current date. • To review the results on a separate PC, for example, in a spreadsheet program, select the ASCII option. • To review the results on the control unit, select the Binary option. 5 Choose the OK button. f After the backup, all measurement results are deleted from the hard disk of the control unit. 6 To remove the USB flash drive, do the following: • Choose the Media Eject button. • Disconnect the USB flash drive from the USB port. When the Reagent Preparing button is purple, yellow, or red, some reagents on board are running short or are empty. You must load the required reagents according to the Reagent Load/Unload List. The color of the Reagent Preparing button indicates the number of available tests or the remaining volume: Reagent below the defined daily requirement Reagent below the defined warning level Reagent empty or missing Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 80 Preparing the analyzer for operation 1 2 r To create a Reagent Load/Unload List 1 During pre-routine operation, select the Preventive Action check box. 1 The purple reagent alarm is only displayed while the Preventive Action check box is selected. 2 Choose the Reagent Preparing button. 3 4 5 3 Choose the Reagent Load/Unload List button and confirm printing. f The report is sent to the print buffer. 4 To view the report, choose the View button. 6 6 In the View button, check which reagents must be unloaded from the analyzer. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Before operation 81 8 7 7 To see the Load List, scroll up. The Load List shows the number of available tests or the remaining reagent volume in mL. 1 The number in brackets indicates the level that is set as daily requirement (purple alarm). 3 Unload and load the reagents according to the report. u Related topics • Loading reagent packs (101) • Unloading reagent packs (104) • Replacing ISE reagents (106) • Changing sample probe detergents (112) • Changing cell detergents (113) Performing recommended calibration and QC Performing recommended calibration and QC 1 Request calibration and QC 2 Print load lists 3 Measure calibrations and QC 4 Check results The analyzer software automatically recommends calibration for all tests that require calibration. Typically, you also perform QC measurements for all tests before operation. In that case you probably do not need to perform a QC measurement during operation. • To activate the recommendation, you must request the calibration and the QC measurement. • You can print load lists, which indicate the required calibrators and controls and their positions on the sample disk. • To start the measurement, you load the calibrators and controls onto the sample disk. • You make sure that all calibrations are valid and all QC results are within the reference range. In this section Requesting recommended calibration and QC (82) Printing load lists for calibrators and controls (83) Measuring calibrators and controls (84) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 82 Preparing the analyzer for operation Requesting recommended calibration and QC QC after calibration The analyzer software recommends calibrations according to the intervals set in the application. However, you must activate the recommended calibrations. To ensure measurement accuracy, you must also perform a QC measurement for all tests. The yellow Calibration And QC Select button indicates that calibration and/or QC are recommended. QC after calibration • Before measuring patient samples j m The analyzer is in Standby status. r To request recommended calibration and QC 1 Choose the Calibration And QC Select button. 1 2 2 In the Calibration group box, choose the Recommended button. f The recommended calibrations are now requested on an internal worldkit and will be performed at the beginning of the next analysis run. 3 Check whether the following steps are necessary: • If a start-up group is defined for this work shift, choose the Start Up button to request the calibration for this group in addition to the recommended tests. • If your work shift starts within less than 24 hours since the last work shift, request an ISE calibration manually. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Before operation 83 5 4 4 In the QC group box, choose the Routine QC button. 5 Choose the Stand By Bottle QC button (optional). f Control measurements are now requested for all reagents on board (in use and standby reagents). u Related topics • About recommended calibration (124) • Performing an ISE calibration (109) Printing load lists for calibrators and controls To see the required calibrators and controls, you can print a Calibrator Load List and QC Load List report. The load lists also show positions on the sample disk assigned to specific calibrators or controls. 1 r To print load lists for calibrators and controls 1 In the Load List group box, select the Calibration option. 2 Choose the Print button. f The Calibration Load List report is printed on paper and sent to the print buffer. It shows all calibrators required on the sample disk. 5 Choose the Close button. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 84 Preparing the analyzer for operation 6 6 On the print-out, check the required calibrators and controls with their assigned sample disk positions.

[illegible]

141 Performing QC after calibration	144 Performing QC for individual standby reagents	146 Checking QC results
149 About the handling of QC results		
150 Workflows for handling and deleting QC results		
152 Transferring QC results		
154 Accumulating QC results		
155 Deleting QC results - without LIS		
156 Deleting QC results - with LIS		
157 Assigning non-barcoded controls to disk positions		
159 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
136 Table of contents Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
3.1 Quality control 137 About quality control (QC) The purpose of QC measurements is to monitor the analyzer performance. QC measurements must be performed in the following cases: • Each day before operation and at the end of operation to ensure measurement accuracy during the entire operation.		
When a new reagent pack is loaded to check the reagent.		
After troubleshooting to check the analyzer performance. In this section About QC workflows (137) About QC types (139) About QC workflows QC measurements during work shift QC measurement is usually performed before operation, before sample processing begins and at the end of operation. However, it can also be done anytime during routine operation. There are different ways to perform QC measurements during a work shift. The workflow shows an example of when you can perform which QC measurement. Before operation QC measurements for all tests During operation QC after calibration At the end of operation QC measurements for individual tests Performing QC measurements for all tests Typically you perform QC measurements for all tests before operation. QC measurements for all tests comprise Routine QC and Standby Bottle QC . If you perform QC measurements for all tests before operation and the QC results are valid, you normally do not need to perform a QC during operation.		
Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
3.1.138 About quality control (QC) Performing QC measurements for individual tests Performing QC after calibration Performing QC measurements for individual standby reagents There are different reasons to perform QC measurements for individual tests: • If a reagent pack is out of range, you must perform QC measurements again. • If you are performing QC after calibration and you want to perform QC measurements for more tests. • If you work with only some tests during operation, you can perform QC measurements for these tests at the end of operation. If you must calibrate a test, you can perform a QC measurement after calibration when suitable controls are loaded on the sample disk. There are different causes to perform QC after calibration: • Before operation, if a calibration is necessary. • During operation, when a new reagent pack with no valid calibration is loaded onto the analyzer. If you are performing QC measurements for individual tests before operation, they have a valid QC. In this case, you do not need to perform a QC measurement during operation when the status of the standby reagent has been changed. There are different causes to perform QC measurements for individual standby reagents: • If a QC result is out of range, you can perform this QC measurement again. • If you are performing QC after calibration and you want to perform QC measurements for more tests including standby reagents. • During operation, when a new reagent pack is loaded as a standby reagent onto the analyzer. u Related topics • About QC types (139) • Performing recommended calibration and QC (81) • Performing QC for individual tests (141) • Performing QC after calibration (144) • Performing QC for individual standby reagents (146) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
3.1 Quality control 139 About QC types Routine QC Standby Bottle QC Timeout QC QC after calibration Depending on the QC type, you can request QC measurements for different tests or for all tests. Each test has one or more controls assigned to it. For QC requests, five QC types are available. With Routine QC, you can request QC measurements for all reagent packs that are in use. Routine QC comprises all activated tests of all installed controls. With Standby Bottle QC, you can request QC measurements for standby reagents. Standby reagents are reagent packs that are already on board but not in use at present. For standby reagents, we recommend requesting QC measurements together with the Routine QC at the beginning of a work shift. QC measurements can be performed at predefined time intervals. When the time interval has elapsed in the System Overview menu the Calib. and QC Select button is displayed in yellow. This indicates a recommended QC measurement. On the QC > Status submenu, the word Timeout in the Cause column also indicates a recommendation for a QC measurement. If the tests are marked with a green bar (Selection column) and the suitable controls are on the sample disk, the analyzer automatically performs necessary QC measurements.		
The controls are measured with the active reagent packs. If the tests are not marked with a green bar, you have to select the tests. For this QC type, no special settings are needed. QC measurements are performed for all tests that are newly calibrated if the suitable controls are on board. You do not have to request this QC. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
3.1.140 About quality control (QC) Manual QC QC measurements can be requested for a single test or multiple tests selected on the QC > Status submenu. For example, if a QC measurement fails, select this control manually on the QC > Status submenu and start a new QC measurement. A manual QC request is indicated by the word Manual in the Cause column on QC > Status submenu. u Related topics • Performing recommended calibration and QC (81) • Performing QC for individual tests (141) • Performing QC after calibration (144) • Performing QC for individual standby reagents (146) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
3.1 Quality control 141 Performing QC for individual tests You can perform QC measurements for individual tests. For example, if you have been working with only some tests during operation, perform QC measurements for these tests at the end of operation.		
This ensures the measurement accuracy for these tests during the entire operation. c Before operation, during operation, or at the end of operation: • If you want to perform QC only for some individual tests. • If a QC result is out of range and you want to perform this QC measurement again. • If you are performing QC after calibration and you want to perform QC measurements for more tests. d m Controls j m The analyzer is in Standby status. r To perform QC measurements for individual tests 1 Choose QC > Status. 1 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
2 Select the desired tests and controls and choose the Select button. • If you want to select all reagent packs that are in use, choose Routine QC Assign. f A green bar is displayed in the Selection column. f Manual is displayed in the Cause column. f The Select button toggles to Deselect. 3 To request the selected controls for measurement, choose the Save button.		
4 Print the QC Load List report. 5 5 WARNING! Incorrect results due to expired or evaporated calibrators or controls. The concentration increases due to evaporation. Only use fresh calibrators and controls. 5 Load the controls onto the sample disk as directed by the QC Load List report.		
Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
3.1.6 Start the measurement. Quality control 143 7 7 Check the results of the QC measurement. f If one QC result is out of range, you must check the control, the analyzer, or the reagent. f If most of the QC results are out of range, you must hold back the sample until the cause for the wrong QC results is known. u Related topics • Printing load lists for calibrators and controls (83) • Measuring calibrators and controls (84) • Checking QC results (149) • List of sample containers (365) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
144 Performing QC after calibration If you must calibrate tests, you can perform QC measurements directly after calibrating these tests.		
These tests, you do not need to request further QC measurements. You may load only suitable controls onto the sample disk. If you must perform QC measurements for more tests, you can perform QC measurements for individual tests 1 Choose QC > Status. 1 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
2 Select the desired tests and controls and choose the Select button. • If you want to select all reagent packs that are in use, choose Routine QC Assign. f A green bar is displayed in the Selection column. f Manual is displayed in the Cause column. f The Select button toggles to Deselect. 3 To request the selected controls for measurement, choose the Save button.		
4 Print the QC Load List report. 5 5 WARNING! Incorrect results due to expired or evaporated calibrators or controls. The concentration increases due to evaporation. Only use fresh calibrators and controls. 5 Load the controls onto the sample disk as directed by the QC Load List report. 5 5 Start the measurement. 6 6 Check the calibration and the results of the QC measurement. u Related topics • Printing load lists for calibrators and controls (83) • Measuring calibrators and controls (84) • Checking calibration results (132) • Checking QC results (149) • List of sample containers (365) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
146 Performing QC for individual standby reagents Performing QC for individual standby reagents You can perform QC measurements for individual standby reagents. For example, if you have loaded a new reagent as a standby reagent onto the analyzer, you must check the reagent with a QC measurement. c Before operation or at the end of operation: • If you want to perform QC only for some individual tests. • If a QC result is out of range, you can perform this QC measurement again. • If you are performing QC after calibration and you want to perform QC measurements for more tests including standby reagents.		
During operation: • If you have loaded a new reagent as a standby reagent onto the analyzer. d m Controls j m The analyzer is in Standby status.		
To perform QC measurements for standby reagents 1 On the QC > Status submenu, choose the Stand By Bottle QC button. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
3 To select all highlighted tests and controls, choose the Select button. f A green bar is displayed in the Selection column. f The Select button toggles to Deselect button. 4 To request the selected controls for measurement, choose the OK button. 5 Print the QC Load List report. 6 6 WARNING! Incorrect results due to expired or evaporated calibrators or controls. The concentration increases due to evaporation. Only use fresh calibrators and controls. 6 Load the controls onto the sample disk as directed by the QC Load List. 7 7 Start the measurement. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
148 Performing QC for individual standby reagents 8 8 Check the results of the QC measurement.		
You can perform this QC measurement again. u Related topics • Measuring calibrators and controls (84) • Checking QC results (149) • List of sample containers (365) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
3.1 Quality control 149 Checking QC results On the QC > Run Status submenu, you can check in an instant whether your QC results are within the corresponding limits. Alternatively, you can view the QC results on the Workplace > Data Review submenu and on the QC > Individual submenu. If the QC results are within the ±2 SD limits, the control, the reagent, and the analyzer function properly. If the QC results fall out the ±2 SD limits, you must check the control, the reagent, and/or the analyzer. 2 3 1 r To check QC results on QC > Run Status 1 Choose QC > Run Status. 2 From the Test drop-down list, choose one of the following: • To view the results of the ±2 500 QC measurements, select the All option. • To view the results of a test, choose the corresponding test. 3 Check if all QC results on a white background (today's data) are within the appropriate SD range (all chart symbols are green). 4 If a QC result falls out of the reference range, perform the following: • Perform the remedies described. • Repeat the QC measurement. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
3.1.150 About the handling of QC results About the handling of QC results The handling of QC results starts after the measurement of QC samples and ends with the deletion of QC results. After the measurement of QC samples, the software displays QC results in different views. You can accumulate the QC results to perform QC statistics on the analyzer. At the end of handling QC results, you must delete the QC results regularly to keep the database compact. Results are displayed in: QC measurements performed Workplace > Data Review (Routine view) QC > Individual Delete Record Transfers one or more QC data to Data Review (QC view) Delete All Transfers all QC data to Data Review (QC view) Workplace > Data Review (QC view) Accumulate Accumulates and transfers one or more QC data to QC > Cumulative QC > Cumulative Delete Record Deletes one or more QC in: • Data Review (QC View) • QC > Individual Delete All Deletes all QC in: • Data Review (QC View) • QC > Individual Delete Deletes one or more QC in: • Data Review (QC View) • QC > Individual Delete Deletes one or more QC in: • QC > Cumulative Displaying of QC results The workflow shows the different submenus that display QC results, and the particular action that moves (deletes or accumulates) the results from one submenu to another submenu.		
When the analyzer has finished the QC measurement, you can view the QC results on the Workplace > Data Review submenu and on the QC > Individual submenu. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual		
3.1 Quality control 151 On the Workplace > Data Review submenu, you can select a control and view the results of the assigned tests. On the QC > Individual submenu, you can select a test and view the result of the assigned control. QC statistics Deletion of QC results If you want to perform the QC statistics on the analyzer, you must accumulate the QC results on the QC > Individual submenu at a regular time interval, for example, at the end of daily operation.		
After accumulation, the corresponding QC results are deleted from the QC > Individual submenu. A new standard deviation (SD) and mean value is calculated and displayed on the QC > Cumulative submenu.		
To keep the performance of the analyzer stable and to minimize the risk of data loss, you must keep the database size compact. Therefore you must delete the QC results regularly on the corresponding submenus. On the QC > Individual submenu and on the Workplace > Data Review submenu, a maximum of 2500 individual QC results can be stored. On the QC > Cumulative submenu, a maximum of 500 cumulated QC results can be stored. Roche Diagnostics cobas c 311 analyzer - Software version		

states. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 287	1 r To clean the detergent aspiration filters 1 Open the left front door. 2 2 Take out the two detergent bottles. 3 3 Remove the tubing from each detergent bottle. 1 CellCln 1 has a filter with a square end. CellCln 2 has a filter with a round end.
5 6 4 Unscrew the filter from each tube end. 5 Wash the filters with deionized water. 6 Then rinse the filters with deionized water. 7 Screw the filters on each tube end. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 288 Monthly maintenance 9	8 Insert the tubing back into each bottle: • Insert each aspiration tube so that the end of the tube touches the bottom of the bottle. • Do not bend the aspiration tubes. I Incorrect results due to incorrect insertion of aspiration tube. If the aspiration tubes are not inserted correctly, the detergent may not be dispensed properly. This may lead to incorrect results. 9 Close the left front door. 10 Perform a cell detergent prime. u Related topics • Performing a cell detergent prime (115) Cleaning the radiator filter To prevent overheating of the analyzer, you must clean the radiator filter of the cooling unit. n Approximately 5 minutes d m Vacuum cleaner m Paper towels m Water for rinsing j m The analyzer is in Standby status.
Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 289	1 r To clean the radiator filter 1 Open the front doors. 2 Vacuum the radiator filter. 3 3 If there is much particulate matter (i.e. dust), clean the filter as follows: • Remove the filter by pulling it forward from the retaining brackets. • Rinse the filter with tap water. • Dry the filter completely with paper towels. • Reinstall the filter. 4 4 Close the front doors. Washing the flow path of concentrated waste solution To prevent clogging of the concentrated waste solution in the flow path, you must wash the flow path. The flow path, the concentrated waste reservoir, and the drainage tubing are cleaned by the following procedure: 1 Choose Utility > Maintenance. 2 Select the (12) Flow Path Wash option. 3 Choose the Select button. 3 2 4 4 Choose the Execute button. Cleaning the water tank To prevent contamination of the entire flow path, you must clean the water tank. This task contains the following procedures: u To disconnect the water tank p (291) u To clean the water tank p (292) u To reconnect the water tank p (293) n Approximately 12 minutes Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 291
d m Lint-free gauze pads m Paper towels m Waste solution receptacle m 0.5% sodium hypochlorite solution m Deionized water j m The analyzer is in Shutdown status. u Shutting down the analyzer (227) 1 r To disconnect the water tank 1 WARNING! Electrical shock by electronic equipment. Cleaning the water tank without turning off the power allows the deionized water supply unit to supply water during cleaning. 1 Shut down the analyzer. Do not reach into the analyzer without switching off power first. 2 2 Turn off the external water supply. 3 3 Open the left front door. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 292 Monthly maintenance 5	4 Close the tap at the outlet of the water tank. 5 To absorb spilled water, place a waste solution receptacle (or paper towels) under the hose unit. 6 6 Separate the joint that connects the water tank and the water hose. 7 7 Remove the float assembly from the tank while pulling out the tank and place the float assembly on a paper towel. 1 r To clean the water tank 1 Visually check the water in the tank. 2 If the water is clear, clean the water tank as follows: • Empty the water from the tank and rinse it with deionized water three times. • Fill up at least 1/3 of the tank of deionized water. • Reconnect the water tank. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 2 Periodic maintenance 293
3 If the water is not clear, perform a thorough cleaning. • Empty the water from the tank and rinse it thoroughly with 0.5% sodium hypochlorite solution. • Use a brush to clean the interior surface. • Wash with deionized water to eliminate the sodium hypochlorite solution. • Fill up at least 1/3 of the tank of deionized water. • Reconnect the water tank. • To reconnect the water tank 1 Wipe the float assembly with gauze pads moistened with deionized water. 2 Return the float assembly and put the tank back in place. 3 3 Reconnect the water hose to the water tank and open the tap. Ensure that all joints are connected properly. 1 If the joints of the water tank are not connected properly, water may leak and damage circuit boards. 4 4 Close the left front door. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 294 Monthly maintenance 5	5 Turn on the external water supply. 6 6 Ensure that the top cover is locked before you switch on the analyzer again. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 295
Every two months maintenance Perform this maintenance tasks at least every two months. In this section Replacing ISE measuring electrodes (Cl, K, Na) (295) Replacing ISE measuring electrodes (Cl, K, Na) (296) u To ensure the measurement accuracy of the ISE unit, you must replace the ISE measuring electrodes. The electrical response level and the slope value (sensitivity) of each measuring electrode slightly decrease with time and use. This task contains the following procedures: u To replace the ISE measuring electrodes p (296) u To perform an ISE check p (299) c You must replace an ISE measuring electrode in the following cases: • The electrode has been in use more than two months. • The test count has reached 9000 tests. • The slope value of the electrode falls outside of the normal range. This is indicated by the data alarms Prep.E or Slop.E. n • Operator time: approximately 10 minutes • System time: approximately 56 minutes (10 minutes for (2) Reagent Prime (IS+REF), 2x 18 minutes for (2) ISE Check, plus 10 minutes in between) d m Sodium (Na) d m Electrode m Potassium (K) d m Electrode m Lint-free gauze pads m Tweezers m Reference solution ISE REF (approx. 1600 µL) m Internal standard ISE IS (approx. 9230 µL) m Calibrators ISE Low (S1), ISE High (S2), and ISE High (S3) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 296 Every two months maintenance j m The analyzer is in Maintenance mode. u Switching the analyzer to Maintenance mode (236) 1 r To replace the ISE measuring electrodes 1 To remove the sipper cover and the cover of the ISE measuring compartment, loosen the screws. 2 2 Disconnect the three electrode wires (Cl, K, and Na). 3 3 To loosen the electrodes from the mounting block, pull the release lever toward the Release position. 1 Injury to your fingers. The release lever is equipped with a spring. Take care of your fingers. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 297	
4 Remove the electrodes from the mounting block using tweezers. 5 5 To ensure measurement precision, thoroughly wipe up any spilled liquid or liquid adhering to connecting parts. 6 6 If an O-ring from the electrodes remains inside the ISE measuring compartment, use tweezers to remove it. 7 7 Remove the black rubber protectors from the electrodes. A Rubber protectors A 8 8 Verify that the connecting part and the new electrodes are provided with O-rings. A A O-ring Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 298 Every two months maintenance 9	9 Insert the new electrodes in the right color-coded position into the ISE measuring compartment while holding the lever in the Release position. 10 10 To fasten the electrodes, set the release lever into the Lock position. 11 11 Reconnect the four color-coded electrode cords with their respective electrodes. 12 12 WARNING! Missing covers of the ISE unit may lead to incorrect results. Always reinstall the covers of the ISE unit after performing maintenance. 12 Reattach the sipper cover and the cover of the ISE measuring compartment. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 299
13 13 Close the top cover of the analyzer and lock it. 14 Terminate the Maintenance mode. 15 Perform a reagent prime with the IS+REF option. 2 3 1 r To perform an ISE check 1 Choose Utility > Maintenance. 2 Select the (2) Check option. 3 Select the (2) ISE Check option. 4 Choose the Select button. 4 6 5 5 Enter "10" cycles. 6 Choose the Execute button. f The electromotive force (EMF) values of the internal standard solution are printed ten times for each electrode. At this point, the results can be ignored. 7 Wait 10 minutes and then repeat the steps 3 to 6. 8 8 Check if any EMF values of the internal standard solution are abnormal: • The difference in successive values for the same electrode should be stable within ± 0.2 mV. • The maximum deviation for Ref EMF over all values should be no more than ± 2 mV. f If an EMF value is abnormal (level error or noise error for example), the corresponding alarm name is printed.	
9 If an EMF value is abnormal, retry ISE check. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 300 Every two months maintenance 10	Perform an ISE calibration before you resume routine operation. • Related topics • Terminating the Maintenance mode (237) • Performing a reagent prime (108) • Performing an ISE calibration (109) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 301
Every three months maintenance Perform this maintenance tasks at least every three months. In this section Cleaning the ultrasonic mixer (301) Replacing the ISE pinch valve tubing and sipper tubing (305) Cleaning the ultrasonic mixer To prevent inadequate mixing and inaccurate results, you must clean the surfaces of the ultrasonic mixer. If you have already cleaned the ultrasonic mixer during monthly maintenance, you do not need to perform this maintenance task separately. u Replacing reaction cells and cleaning incubator bath (278) This task contains the following procedures: u To clean the surfaces of the ultrasonic mixer p (302) u To check the intensity of the ultrasonic output p (305) n Approximately 15 minutes d m Cotton swabs m Paper towels m EcoTergent reagent pack (surfactant for the incubator, consumption 4.3 mL) m Deionized water j m The analyzer is in Standby status. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 302 Every three months maintenance 3	2 1 r To clean the surfaces of the ultrasonic mixer 1 Choose Utility > Maintenance. 2 Select the (9) Incubator Bath Cleaning option. 3 Choose the Select button. 4 4 Choose the Execute button. f The water is drained from the incubator bath and the control unit is powered off (shutdown). 5 5 After the computer power supply has been powered off, switch off the power of the analytical unit. 1 Do not open the top cover until the analyzer is in shutdown. 6 6 To remove the cell cover attached to the ultrasonic mixer, loosen the screws. A A Screw Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 303
7 7 To remove the segment of the reaction cells near the ultrasonic mixer, remove the respective screws. 1 Do not touch the surfaces of the reaction cells. 8 8 Wipe the polished surfaces of the ultrasonic mixer with a cotton swab moistened with deionized water. A Ultrasonic mixer A 9 9 Put the removed segment of reaction cells back in its previous place. 10 10 Reinstall the cell cover. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 304 Every three months maintenance 11	11 Close the top cover of the analyzer and lock it. 12 12 Switch on the analyzer. f The incubator bath is filled with some water. A yellow alarm is issued, indicating that the level of incubation water is below the lower limit. 14 13 13 Choose Utility > Maintenance. 14 Select the (4) Incubation Water Exchange option. 15 Choose the Select button. 15 16 16 Choose the Execute button. f The incubator is exchanged and EcoTergent (4.3 mL) is added. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 305
2 3 1 r To check the intensity of the ultrasonic output 1 Choose Utility > Maintenance. 2 Select the (2) Check option. 3 Select the (16) Cuvette Mixing option. 4 Choose the Select button. 4 6 5 5 Enter "10" cycles. 6 Choose the Execute button. 7 If an alarm occurs, follow the remedies described in the alarm message. u Related topics • Replacing reaction cells and cleaning incubator bath (278) Replacing the ISE pinch valve tubing and sipper tubing After long use, the ISE pinch valve tubing and the sipper tubing will gradually wear out. To ensure the accuracy of sample aspiration, you must replace both tubings. This task contains the following procedures: u To replace the ISE pinch valve tubing and the sipper tubing p (306) u To perform an ISE check p (308) n • Operator time: approximately 9 minutes • System time: approximately 46 minutes (2x 18 minutes for (2) ISE Check, plus 10 minutes in between) d m ISE pinch valve tubing m ISE sipper tubing m Calibrators ISE Low (S1), ISE High (S2), and ISE High (S3) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 306 Every three months maintenance j m The analyzer is in Maintenance mode. u Switching the analyzer to Maintenance mode (236) 1 r To replace the ISE pinch valve tubing and the sipper tubing 1 To remove the sipper cover and the cover of the ISE measuring compartment, loosen the screws. 2 2 Carefully remove both ends of the ISE pinch valve tubing from their tube connectors and pull the tubing through the pinch valve. A ISE pinch valve tubing A 3 3 Insert the new pinch valve tubing through the pinch valve and attach both ends of the tubing to the connectors. 1 Make sure that there is no slack in the tubing. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 307	
4 4 Carefully remove both ends of the ISE sipper tubing from the tube connectors. A ISE sipper tubing A 5 5 Attach both ends of the new sipper tubing to the tube connectors. 6 6 WARNING! Missing covers of the ISE unit may lead to incorrect results. Always reinstall the covers of the ISE unit after performing maintenance. 6 Reattach the sipper cover and the cover of the ISE measuring compartment. 7 7 Close the top cover of the analyzer and lock it. 8 Terminate the Maintenance mode. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 308 Every three months maintenance 2	3 1 r To perform an ISE check 1 Choose Utility > Maintenance. 2 Select the (2) Check option. 3 Select the (2) ISE Check option. 4 Choose the Select button. 4 6 5 5 Enter "10" cycles. 6 Choose the Execute button. f The electromotive force (EMF) values of the internal standard solution are printed ten times for each electrode. At this point, the results can be ignored. 7 Wait 10 minutes and then repeat the steps 3 to 6. 8 8 Check if any EMF values of the internal standard solution are abnormal: • The difference in successive values for the same electrode should be stable within ± 0.2 mV. • Normal IS EMF values lie between NA: -10 to -90 mV, K -10 to -90 mV and Cl 80-160 mV, Ref -7 to 7 mV. • The maximum deviation for Ref EMF over all values should be no more than ± 2 mV. f If an EMF value is abnormal (level error or noise error for example), the corresponding alarm name is printed. 9 If an EMF value is abnormal, retry ISE check. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10
Perform an ISE calibration before you resume routine operation. u Related topics • Terminating the Maintenance mode (237) • Performing an ISE calibration (109) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 309	Every six months maintenance Perform this maintenance tasks at least every six months. In this section Cleaning the inlet water filter (309) Cleaning the cooling fan (312) Replacing the photometer lamp (313) Replacing the ISE reference electrode (320) Cleaning the inlet water filter To prevent clogging of the water system, you must clean the inlet water filter. n Approximately 5 minutes d m Water pump filter (in case replacement is required) m Paper towels m Beaker or container, 500 mL m Deionized water j m The analyzer is in Shutdown status. u Shutting down the analyzer (227) 1 To clean the inlet water filter 1 Turn off the external water supply. 1 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 310 Every six months maintenance 3
2 2 Place a beaker (500 mL) beneath the inlet water manifold. 3 3 Turn the ring on the water filter cap counterclockwise and disconnect the inlet water hose. 4 4 Remove the water filter and insert the hose into the container. 5 5 Clean the filter thoroughly with deionized water. 6 6 Reinstall the filter. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 311	7 7 CAUTION! Incorrect results due to loose water filter cap. If the inlet water hose is not reconnected to the inlet water manifold correctly, it may cause instrument malfunction. This can lead to incorrect measurement results. 7 Reconnect the inlet water hose to the inlet water manifold. 8 8 Turn on the external water supply. 9 9 Ensure that the water filter is locked on the analyzer again. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 312 Every six months maintenance Cleaning the cooling fan To remove dust and dirt, you must clean the cooling fan at the rear of the analyzer. 1 CAUTION Personal injury through touching the fan There is a grill to protect the operator from coming into contact with the moving fan. If the fan is cleaned when the analyzer is not in Shutdown status, there is a danger of injury. r Shut down the analyzer before cleaning. r Do not insert your fingers in the openings of the ventilation grill. n Approximately 5 minutes d m Vacuum cleaner j m The analyzer is in Shutdown status. u Shutting down the analyzer (227) 1 r To clean the cooling fan 1 Vacuum dust, dirt, and other debris from the cooling fan at the back of the analyzer. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 313
Replacing the photometer lamp To ensure the reproducibility of measurement, you must replace the photometer lamp. After installing a new photometer lamp, you must wait 30 minutes for the lamp to stabilize. Then you must perform a cell blank measurement. This is necessary to compensate for a potential change in light intensity. q We recommend combining this maintenance task with the monthly cleaning of the incubator bath and with the cleaning of the ultrasonic mixer every three months. You must check the light intensity of the photometer lamp daily to maintain reproducibility of measurement. The photometer check is included in the daily wake-up pipe. u Checking the analyzer and performing maintenance (66) u List of recommended maintenance pipes (244) f If the photometer check value exceeds 14 000 absorbance units at any wavelength, check the following points before replacing the photometer lamp: • Make sure that the reaction cells, through the pinch valve, A ISE pinch valve tubing A 1 A ISE pinch valve tubing and sipper tubing After long use, the ISE pinch valve tubing and the sipper tubing will gradually wear out. To ensure the accuracy of sample aspiration, you must replace both tubings. This task contains the following procedures: u To replace the ISE pinch valve tubing and the sipper tubing p (306) u To perform an ISE check p (308) n • Operator time: approximately 9 minutes • System time: approximately 46 minutes (2x 18 minutes for (2) ISE Check, plus 10 minutes in between) d m ISE pinch valve tubing m ISE sipper tubing m Calibrators ISE Low (S1), ISE High (S2), and ISE High (S3) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 306 Every three months maintenance j m The analyzer is in Maintenance mode. u Switching the analyzer to Maintenance mode (236) 1 r To replace the ISE pinch valve tubing and the sipper tubing 1 To remove the sipper cover and the cover of the ISE measuring compartment, loosen the screws. 2 2 Carefully remove both ends of the ISE pinch valve tubing from their tube connectors and pull the tubing through the pinch valve. A ISE pinch valve tubing A 3 3 Insert the new pinch valve tubing through the pinch valve and attach both ends of the tubing to the connectors. 1 Make sure that there is no slack in the tubing. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 307	
4 4 Carefully remove both ends of the ISE sipper tubing from the tube connectors. A ISE sipper tubing A 5 5 Attach both ends of the new sipper tubing to the tube connectors. 6 6 WARNING! Missing covers of the ISE unit may lead to incorrect results. Always reinstall the covers of the ISE unit after performing maintenance. 6 Reattach the sipper cover and the cover of the ISE measuring compartment. 7 7 Close the top cover of the analyzer and lock it. 8 Terminate the Maintenance mode. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 308 Every three months maintenance 2	3 1 r To perform an ISE check 1 Choose Utility > Maintenance. 2 Select the (2) Check option. 3 Select the (2) ISE Check option. 4 Choose the Select button. 4 6 5 5 Enter "10" cycles. 6 Choose the Execute button. f The electromotive force (EMF) values of the internal standard solution are printed ten times for each electrode. At this point, the results can be ignored. 7 Wait 10 minutes and then repeat the steps 3 to 6. 8 8 Check if any EMF values of the internal standard solution are abnormal: • The difference in successive values for the same electrode should be stable within ± 0.2 mV. • Normal IS EMF values lie between NA: -10 to -90 mV, K -10 to -90 mV and Cl 80-160 mV, Ref -7 to 7 mV. • The maximum deviation for Ref EMF over all values should be no more than ± 2 mV. f If an EMF value is abnormal (level error or noise error for example), the corresponding alarm name is printed. 9 If an EMF value is abnormal, retry ISE check. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10
Perform an ISE calibration before you resume routine operation. u Related topics • Terminating the Maintenance mode (237) • Performing a reagent prime (108) • Performing an ISE calibration (109) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10	Perform an ISE calibration before you resume routine operation. u Related topics • Terminating the Maintenance mode (237) • Performing a reagent prime (108) • Performing an ISE calibration (109) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10
Every six months maintenance Perform this maintenance tasks at least every six months. In this section Cleaning the inlet water filter (309) Cleaning the cooling fan (312) Replacing the photometer lamp (313) Replacing the ISE reference electrode (320) Cleaning the inlet water filter To prevent clogging of the water system, you must clean the inlet water filter. n Approximately 5 minutes d m Water pump filter (in case replacement is required) m Paper towels m Beaker or container, 500 mL m Deionized water j m The analyzer is in Shutdown status. u Shutting down the analyzer (227) 1 To clean the inlet water filter 1 Turn off the external water supply. 1 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 310 Every six months maintenance 3	2 2 Place a beaker (500 mL) beneath the inlet water manifold. 3 3 Turn the ring on the water filter cap counterclockwise and disconnect the inlet water hose. 4 4 Remove the water filter and insert the hose into the container. 5 5 Clean the filter thoroughly with deionized water. 6 6 Reinstall the filter. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 311
7 7 CAUTION! Incorrect results due to loose water filter cap. If the inlet water hose is not reconnected to the inlet water manifold correctly, it may cause instrument malfunction. This can lead to incorrect measurement results. 7 Reconnect the inlet water hose to the inlet water manifold. 8 8 Turn on the external water supply. 9 9 Ensure that the water filter is locked on the analyzer again. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 312 Every six months maintenance Cleaning the cooling fan To remove dust and dirt, you must clean the cooling fan at the rear of the analyzer. 1 CAUTION Personal injury through touching the fan There is a grill to protect the operator from coming into contact with the moving fan. If the fan is cleaned when the analyzer is not in Shutdown status, there is a danger of injury. r Shut down the analyzer before cleaning. r Do not insert your fingers in the openings of the ventilation grill. n Approximately 5 minutes d m Vacuum cleaner j m The analyzer is in Shutdown status. u Shutting down the analyzer (227) 1 r To clean the cooling fan 1 Vacuum dust, dirt, and other debris from the cooling fan at the back of the analyzer. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 313	Replacing the photometer lamp To ensure the reproducibility of measurement, you must replace the photometer lamp. After installing a new photometer lamp, you must wait 30 minutes for the lamp to stabilize. Then you must perform a cell blank measurement. This is necessary to compensate for a potential change in light intensity. q We recommend combining this maintenance task with the monthly cleaning of the incubator bath and with the cleaning of the ultrasonic mixer every three months. You must check the light intensity of the photometer lamp daily to maintain reproducibility of measurement. The photometer check is included in the daily wake-up pipe. u Checking the analyzer and performing maintenance (66) u List of recommended maintenance pipes (244) f If the photometer check value exceeds 14 000 absorbance units at any wavelength, check the following points before replacing the photometer lamp: • Make sure that the reaction cells, through the pinch valve, A ISE pinch valve tubing A 1 A ISE pinch valve tubing and sipper tubing After long use, the ISE pinch valve tubing and the sipper tubing will gradually wear out. To ensure the accuracy of sample aspiration, you must replace both tubings. This task contains the following procedures: u To replace the ISE pinch valve tubing and the sipper tubing p (306) u To perform an ISE check p (308) n • Operator time: approximately 9 minutes • System time: approximately 46 minutes (2x 18 minutes for (2) ISE Check, plus 10 minutes in between) d m ISE pinch valve tubing m ISE sipper tubing m Calibrators ISE Low (S1), ISE High (S2), and ISE High (S3) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 306 Every three months maintenance j m The analyzer is in Maintenance mode. u Switching the analyzer to Maintenance mode (236) 1 r To replace the ISE pinch valve tubing and the sipper tubing 1 To remove the sipper cover and the cover of the ISE measuring compartment, loosen the screws. 2 2 Carefully remove both ends of the ISE pinch valve tubing from their tube connectors and pull the tubing through the pinch valve. A ISE pinch valve tubing A 3 3 Insert the new pinch valve tubing through the pinch valve and attach both ends of the tubing to the connectors. 1 Make sure that there is no slack in the tubing. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 307
4 4 Carefully remove both ends of the ISE sipper tubing from the tube connectors. A ISE sipper tubing A 5 5 Attach both ends of the new sipper tubing to the tube connectors. 6 6 WARNING! Missing covers of the ISE unit may lead to incorrect results. Always reinstall the covers of the ISE unit after performing maintenance. 6 Reattach the sipper cover and the cover of the ISE measuring compartment. 7 7 Close the top cover of the analyzer and lock it. 8 Terminate the Maintenance mode. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 308 Every three months maintenance 2	3 1 r To perform an ISE check 1 Choose Utility > Maintenance. 2 Select the (2) Check option. 3 Select the (2) ISE Check option. 4 Choose the Select button. 4 6 5 5 Enter "10" cycles. 6 Choose the Execute button. f The electromotive force (EMF) values of the internal standard solution are printed ten times for each electrode. At this point, the results can be ignored. 7 Wait 10 minutes and then repeat the steps 3 to 6. 8 8 Check if any EMF values of the internal standard solution are abnormal: • The difference in successive values for the same electrode should be stable within ± 0.2 mV. • Normal IS EMF values lie between NA: -10 to -90 mV, K -10 to -90 mV and Cl 80-160 mV, Ref -7 to 7 mV. • The maximum deviation for Ref EMF over all values should be no more than ± 2 mV. f If an EMF value is abnormal (level error or noise error for example), the corresponding alarm name is printed. 9 If an EMF value is abnormal, retry ISE check. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10
Perform an ISE calibration before you resume routine operation. u Related topics • Terminating the Maintenance mode (237) • Performing a reagent prime (108) • Performing an ISE calibration (109) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10	Perform an ISE calibration before you resume routine operation. u Related topics • Terminating the Maintenance mode (237) • Performing a reagent prime (108) • Performing an ISE calibration (109) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10
Every six months maintenance Perform this maintenance tasks at least every six months. In this section Cleaning the inlet water filter (309) Cleaning the cooling fan (312) Replacing the photometer lamp (313) Replacing the ISE reference electrode (320) Cleaning the inlet water filter To prevent clogging of the water system, you must clean the inlet water filter. n Approximately 5 minutes d m Water pump filter (in case replacement is required) m Paper towels m Beaker or container, 500 mL m Deionized water j m The analyzer is in Shutdown status. u Shutting down the analyzer (227) 1 To clean the inlet water filter 1 Turn off the external water supply. 1 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 310 Every six months maintenance 3	2 2 Place a beaker (500 mL) beneath the inlet water manifold. 3 3 Turn the ring on the water filter cap counterclockwise and disconnect the inlet water hose. 4 4 Remove the water filter and insert the hose into the container. 5 5 Clean the filter thoroughly with deionized water. 6 6 Reinstall the filter. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 311
7 7 CAUTION! Incorrect results due to loose water filter cap. If the inlet water hose is not reconnected to the inlet water manifold correctly, it may cause instrument malfunction. This can lead to incorrect measurement results. 7 Reconnect the inlet water hose to the inlet water manifold. 8 8 Turn on the external water supply. 9 9 Ensure that the water filter is locked on the analyzer again. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 312 Every six months maintenance Cleaning the cooling fan To remove dust and dirt, you must clean the cooling fan at the rear of the analyzer. 1 CAUTION Personal injury through touching the fan There is a grill to protect the operator from coming into contact with the moving fan. If the fan is cleaned when the analyzer is not in Shutdown status, there is a danger of injury. r Shut down the analyzer before cleaning. r Do not insert your fingers in the openings of the ventilation grill. n Approximately 5 minutes d m Vacuum cleaner j m The analyzer is in Shutdown status. u Shutting down the analyzer (227) 1 r To clean the cooling fan 1 Vacuum dust, dirt, and other debris from the cooling fan at the back of the analyzer. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 313	Replacing the photometer lamp To ensure the reproducibility of measurement, you must replace the photometer lamp. After installing a new photometer lamp, you must wait 30 minutes for the lamp to stabilize. Then you must perform a cell blank measurement. This is necessary to compensate for a potential change in light intensity. q We recommend combining this maintenance task with the monthly cleaning of the incubator bath and with the cleaning of the ultrasonic mixer every three months. You must check the light intensity of the photometer lamp daily to maintain reproducibility of measurement. The photometer check is included in the daily wake-up pipe. u Checking the analyzer and performing maintenance (66) u List of recommended maintenance pipes (244) f If the photometer check value exceeds 14 000 absorbance units at any wavelength, check the following points before replacing the photometer lamp: • Make sure that the reaction cells, through the pinch valve, A ISE pinch valve tubing A 1 A ISE pinch valve tubing and sipper tubing After long use, the ISE pinch valve tubing and the sipper tubing will gradually wear out. To ensure the accuracy of sample aspiration, you must replace both tubings. This task contains the following procedures: u To replace the ISE pinch valve tubing and the sipper tubing p (306) u To perform an ISE check p (308) n • Operator time: approximately 9 minutes • System time: approximately 46 minutes (2x 18 minutes for (2) ISE Check, plus 10 minutes in between) d m ISE pinch valve tubing m ISE sipper tubing m Calibrators ISE Low (S1), ISE High (S2), and ISE High (S3) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 306 Every three months maintenance j m The analyzer is in Maintenance mode. u Switching the analyzer to Maintenance mode (236) 1 r To replace the ISE pinch valve tubing and the sipper tubing 1 To remove the sipper cover and the cover of the ISE measuring compartment, loosen the screws. 2 2 Carefully remove both ends of the ISE pinch valve tubing from their tube connectors and pull the tubing through the pinch valve. A ISE pinch valve tubing A 3 3 Insert the new pinch valve tubing through the pinch valve and attach both ends of the tubing to the connectors. 1 Make sure that there is no slack in the tubing. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 307
4 4 Carefully remove both ends of the ISE sipper tubing from the tube connectors. A ISE sipper tubing A 5 5 Attach both ends of the new sipper tubing to the tube connectors. 6 6 WARNING! Missing covers of the ISE unit may lead to incorrect results. Always reinstall the covers of the ISE unit after performing maintenance. 6 Reattach the sipper cover and the cover of the ISE measuring compartment. 7 7 Close the top cover of the analyzer and lock it. 8 Terminate the Maintenance mode. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 308 Every three months maintenance 2	3 1 r To perform an ISE check 1 Choose Utility > Maintenance. 2 Select the (2) Check option. 3 Select the (2) ISE Check option. 4 Choose the Select button. 4 6 5 5 Enter "10" cycles. 6 Choose the Execute button. f The electromotive force (EMF) values of the internal standard solution are printed ten times for each electrode. At this point, the results can be ignored. 7 Wait 10 minutes and then repeat the steps 3 to 6. 8 8 Check if any EMF values of the internal standard solution are abnormal: • The difference in successive values for the same electrode should be stable within ± 0.2 mV. • Normal IS EMF values lie between NA: -10 to -90 mV, K -10 to -90 mV and Cl 80-160 mV, Ref -7 to 7 mV. • The maximum deviation for Ref EMF over all values should be no more than ± 2 mV. f If an EMF value is abnormal (level error or noise error for example), the corresponding alarm name is printed. 9 If an EMF value is abnormal, retry ISE check. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10
Perform an ISE calibration before you resume routine operation. u Related topics • Terminating the Maintenance mode (237) • Performing a reagent prime (108) • Performing an ISE calibration (109) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10	Perform an ISE calibration before you resume routine operation. u Related topics • Terminating the Maintenance mode (237) • Performing a reagent prime (108) • Performing an ISE calibration (109) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10
Every six months maintenance Perform this maintenance tasks at least every six months. In this section Cleaning the inlet water filter (309) Cleaning the cooling fan (312) Replacing the photometer lamp (313) Replacing the ISE reference electrode (320) Cleaning the inlet water filter To prevent clogging of the water system, you must clean the inlet water filter. n Approximately 5 minutes d m Water pump filter (in case replacement is required) m Paper towels m Beaker or container, 500 mL m Deionized water j m The analyzer is in Shutdown status. u Shutting down the analyzer (227) 1 To clean the inlet water filter 1 Turn off the external water supply. 1 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 310 Every six months maintenance 3	2 2 Place a beaker (500 mL) beneath the inlet water manifold. 3 3 Turn the ring on the water filter cap counterclockwise and disconnect the inlet water hose. 4 4 Remove the water filter and insert the hose into the container. 5 5 Clean the filter thoroughly with deionized water. 6 6 Reinstall the filter. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 311
7 7 CAUTION! Incorrect results due to loose water filter cap. If the inlet water hose is not reconnected to the inlet water manifold correctly, it may cause instrument malfunction. This can lead to incorrect measurement results. 7 Reconnect the inlet water hose to the inlet water manifold. 8 8 Turn on the external water supply. 9 9 Ensure that the water filter is locked on the analyzer again. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 312 Every six months maintenance Cleaning the cooling fan To remove dust and dirt, you must clean the cooling fan at the rear of the analyzer. 1 CAUTION Personal injury through touching the fan There is a grill to protect the operator from coming into contact with the moving fan. If the fan is cleaned when the analyzer is not in Shutdown status, there is a danger of injury. r Shut down the analyzer before cleaning. r Do not insert your fingers in the openings of the ventilation grill. n Approximately 5 minutes d m Vacuum cleaner j m The analyzer is in Shutdown status. u Shutting down the analyzer (227) 1 r To clean the cooling fan 1 Vacuum dust, dirt, and other debris from the cooling fan at the back of the analyzer. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 313	Replacing the photometer lamp To ensure the reproducibility of measurement, you must replace the photometer lamp. After installing a new photometer lamp, you must wait 30 minutes for the lamp to stabilize. Then you must perform a cell blank measurement. This is necessary to compensate for a potential change in light intensity. q We recommend combining this maintenance task with the monthly cleaning of the incubator bath and with the cleaning of the ultrasonic mixer every three months. You must check the light intensity of the photometer lamp daily to maintain reproducibility of measurement. The photometer check is included in the daily wake-up pipe. u Checking the analyzer and performing maintenance (66) u List of recommended maintenance pipes (244) f If the photometer check value exceeds 14 000 absorbance units at any wavelength, check the following points before replacing the photometer lamp: • Make sure that the reaction cells, through the pinch valve, A ISE pinch valve tubing A 1 A ISE pinch valve tubing and sipper tubing After long use, the ISE pinch valve tubing and the sipper tubing will gradually wear out. To ensure the accuracy of sample aspiration, you must replace both tubings. This task contains the following procedures: u To replace the ISE pinch valve tubing and the sipper tubing p (306) u To perform an ISE check p (308) n • Operator time: approximately 9 minutes • System time: approximately 46 minutes (2x 18 minutes for (2) ISE Check, plus 10 minutes in between) d m ISE pinch valve tubing m ISE sipper tubing m Calibrators ISE Low (S1), ISE High (S2), and ISE High (S3) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 306 Every three months maintenance j m The analyzer is in Maintenance mode. u Switching the analyzer to Maintenance mode (236) 1 r To replace the ISE pinch valve tubing and the sipper tubing 1 To remove the sipper cover and the cover of the ISE measuring compartment, loosen the screws. 2 2 Carefully remove both ends of the ISE pinch valve tubing from their tube connectors and pull the tubing through the pinch valve. A ISE pinch valve tubing A 3 3 Insert the new pinch valve tubing through the pinch valve and attach both ends of the tubing to the connectors. 1 Make sure that there is no slack in the tubing. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 307
4 4 Carefully remove both ends of the ISE sipper tubing from the tube connectors. A ISE sipper tubing A 5 5 Attach both ends of the new sipper tubing to the tube connectors. 6 6 WARNING! Missing covers of the ISE unit may lead to incorrect results. Always reinstall the covers of the ISE unit after performing maintenance. 6 Reattach the sipper cover and the cover of the ISE measuring compartment. 7 7 Close the top cover of the analyzer and lock it. 8 Terminate the Maintenance mode. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 308 Every three months maintenance 2	3 1 r To perform an ISE check 1 Choose Utility > Maintenance. 2 Select the (2) Check option. 3 Select the (2) ISE Check option. 4 Choose the Select button. 4 6 5 5 Enter "10" cycles. 6 Choose the Execute button. f The electromotive force (EMF) values of the internal standard solution are printed ten times for each electrode. At this point, the results can be ignored. 7 Wait 10 minutes and then repeat the steps 3 to 6. 8 8 Check if any EMF values of the internal standard solution are abnormal: • The difference in successive values for the same electrode should be stable within ± 0.2 mV. • Normal IS EMF values lie between NA: -10 to -90 mV, K -10 to -90 mV and Cl 80-160 mV, Ref -7 to 7 mV. • The maximum deviation for Ref EMF over all values should be no more than ± 2 mV. f If an EMF value is abnormal (level error or noise error for example), the corresponding alarm name is printed. 9 If an EMF value is abnormal, retry ISE check. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10
Perform an ISE calibration before you resume routine operation. u Related topics • Terminating the Maintenance mode (237) • Performing a reagent prime (108) • Performing an ISE calibration (109) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10	Perform an ISE calibration before you resume routine operation. u Related topics • Terminating the Maintenance mode (237) • Performing a reagent prime (108) • Performing an ISE calibration (109) Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 309 Every six months maintenance 10
Every six months maintenance Perform this maintenance tasks at least every six months. In this section Cleaning the inlet water filter (309) Cleaning the cooling fan (312) Replacing the photometer lamp (313) Replacing the ISE reference electrode (320) Cleaning the inlet water filter To prevent clogging of the water system, you must clean the inlet water filter. n Approximately 5 minutes d m Water pump filter (in case replacement is required) m Paper towels m Beaker or container, 500 mL m Deionized water j m The analyzer is in Shutdown status. u Shutting down the analyzer (227) 1 To clean the inlet water filter 1 Turn off the external water supply. 1 Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 310 Every six months maintenance 3	2 2 Place a beaker (500 mL) beneath the inlet water manifold. 3 3 Turn the ring on the water filter cap counterclockwise and disconnect the inlet water hose. 4 4 Remove the water filter and insert the hose into the container. 5 5 Clean the filter thoroughly with deionized water. 6 6 Reinstall the filter. Roche Diagnostics cobas c 311 analyzer - Software version 01-10 - Operator's Manual - 3.1 Periodic maintenance 311
7	

180 µ L, in 1 µ L increments (5-19 µ L + 20 µ L water) 3 timings possible R1 R2 0 min R3 Up to 42 3.5 min Number of cobas c packs on board cobas storage Remaining reagent volume cobas time 5-15 °C Automatic test countdown with each pipetting y Reagent system of the photometric unit Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Specifications 363 List of ISE unit specifications General characteristics The ISE unit is located next to the reaction disk. It measures electrolytes with an ion-selective electrode. The reagent system is shared with the photometric unit. The characteristics of the reagent system for the ISE unit are given in this topic. The ISE unit has the following general characteristics. Applications Detection system Measuring range (serum) Measuring range (urine) Measuring range (urine rerun with reduced sample volume) Characteristic Cl⁻ Na⁺ + K⁺ + Cl⁻ K⁺ + Cl⁻ Na⁺ + K⁺ + Na⁺ + K⁺ + , Cl⁻ in serum and urine Ion-selective electrode system Na⁺ + 80-180 mmol/L 1.5-10.0 mmol/L 60-140 mmol/L 10-250 mmol/L 1-100 mmol/L 10-250 mmol/L 250-375 mmol/L 100-150 mmol/L 250-375 mmol/L 37 °C ± 2 °C (during operation ± 0.5 °C) 24 s Capacitance sensing technology Measuring temperature Cycle time Liquid level sensor y ISE unit Reagent consumption Bottle volume Reagent system The reagent system has the following characteristics for the ISE unit. Characteristic ISE DIL ISE IS ISE IS ISE DIL ISE IS ISE IS 355 µ L/sample (dilution to 1/31) Measurements in succession Single measurement Measurements in succession Single measurement 600 mL 300 mL 300 mL 710 µ L/sample 1420 µ L/sample 130 µ L/sample 195 µ L/sample y Reagent system of the ISE unit Volume Aspirated from reagent bottle Discharged to reaction cell (1) Aspirated by the ISE sipper ISE IS 710 µ L 450 µ L ISE DIL 355 µ L 291 µ L 400 µ L 250 µ L (2) y Reagent pipetting volumes during analysis (1) The discharged volumes are lower because of the required dummy volume of the reagent pipetter. ISE REF - - 65 µ L Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 364 Analytical specifications (2) Diluted sample. Dilution ratio 1:31 Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Specifications 365 List of sample containers The sample containers hold the undiluted or manually prediluted samples as well as calibrators and controls on the sample disk. Standard cups and micro cups are both called cups . Cups can be inserted into sample tubes with 16 mm diameter, which is called cup on tube . Cups can also be used without tubes. A A Sample tube B Micro cup B C D C Standard cup D False bottom tube Container Primary sample tube Diameter x length 16 mm x 100 mm 16 mm x 75 mm 13 mm x 100 mm 13 mm x 100 mm Dead volume 1000 µ L 1000 µ L 500 µ L 500 µ L 100 µ L Hitachi standard cup or cobas sample cup, 2.5 mL Hitachi micro cup (1) , 1.5 mL Roche Diagnostics Standard False Bottom Tube (FBT) (3) Freely definable type of false bottom tube / non-standard tube (3) 13 mm x 75 mm 12-16 mm x 73-102 mm 50 µ L (2) L y Accepted sample containers (1) Not applicable to calibrators and controls (2) A dead volume of 50 µ L is available if the Hitachi micro cup is directly loaded on the sample disk (not as a cup on tube) and the sample cup size is set to Micro on the Workplace > Test Selection screen. (3) Only one type of false bottom / non-standard tube can be used on the analyzer at the same time. q Avoid overfilling sample tubes and cups.

Leave approx. 10 mm clearance between the liquid and the top of the tube or cup. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 366 List of barcode specifications List of barcode specifications Sample IDs or Patient IDs are encoded in barcodes. Barcodes are placed on sample containers and read by barcode readers next to the inner and outer ring of the sample disk. Barcodes used with the cobas c 311 analyzer must be in compliance with one of the following standards: • NW7 (Codabar) • Code 39 • ITF • Code 128 Sample barcodes The following specifications apply to the various barcode types: Reading method Used barcode symbol Check digit Number of ID digits Characters usable for ID Check digit Specification Scanning with CCD sensor o NW7 (Codabar) o Code 39 o ITF o Code 128 Must be used to prevent scanning errors. NW7 Code 39 ITF Code 128 NW7 Code 39 ITF Code 128 2-22 digits + 1 digit (with check digit) 2-22 digits + 1 digit (with check digit) 3-21 digits + 1 digit (with check digit) 3-22 digits + 2 digit (with check digit) 0 to 9 (without check digit) 0 to 9, -, /, ., , \$, %, + (with check digit) 0 to 9, A to Z, -, ., [, /, +, \$, % 0 to 9 Alphanumerics (excluding those alphanumerics assigned to functions and communications) NW7 Code 39 ITF Code 128 Modulus 16 Modulus 11 Modulus 10/2 weight Modulus 10/3 weight 7 check DR weighted Modulus 11 Modulus 4/3 Modulus 10/3 weight Modulus 103 y Sample barcode specifications Specifications of barcode labels The barcode labels must comply with the following specifications: Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Specifications 367 Secure a margin of 5 mm or more on each edge of barcode label. •PP •PP IRUPPVDPDSOHWXEHV •PP IRUPPVDPDSOHWXEHV •PP •PP •PP Bar and space Reflectance and PCS value Colors Other y Specifications of barcode labels Specification Minimum bar (space) width Ratio of narrow bar (space) to wide bar (space) Space reflectance PCS value Bar Space Bar print must not be blurred. Label indications must not be contaminated. 0.19 up to 0.20 mm 1.2/5 to 3 ≥ 75% ≥ 0.7 Black White u Related topics • Placing barcode labels on sample tubes (213) Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 368 List of accessories and consumables List of accessories and consumables Image of product This is a list of globally available accessories and consumables.

For ordering information, contact your local sales representative. Product name Reaction cells ISE cleaning solution (SysClean) ECO-D solution (EcoTergent) Concentrated liquid waste container y List of consumables u Related topics • About positions of reagents (95) • List of tools, materials, and solutions (254) • List of spare parts and replacement intervals (255) • List of sample containers (365) Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Glossary 369 Alarm button - maintenance mode Glossary This glossary explains technical terms used in this manual. Alarm button Global button used to display the list of current alarms. If an alarm was issued, the button flashes yellow or red, automatic or red, without operator intervention, of tests that have results with data alarms.

auxiliary reagent A non-test-specific reagent for washing and rinsing reaction cells, sample probes, and reagent probes. calibration masking A function that masks a reagent when no valid calibration is available. cobas e-LabPerformance Application that uploads statistical data from the instruments via cobas link to the remote service platform for performance monitoring and QC management. cobas e-library An electronic archive on the cobas link data station containing e-library package inserts and e-barcodes.

For each application, calibrator, and control, the archive contains the latest version of the respective e-library package and the preceding version. The status line is yellow and displays the Standby status. cobas e-services A set of cobas link applications. The set includes: cobas e-library, cobas e-support, and cobas e-LabPerformance. cobas e-support Screen sharing via cobas link with a Roche Service representative for troubleshooting and support. cobas link Software package running on a dedicated computer in the laboratory, the cobas link data station. The data station and the software are together referred to as a cobas link. A cobas link connects one or more cobas instruments to the remote service platform.

cobas link certificate A 16-digit PIN code issued for a specific cobas link. The certificate is valid for 14 days, i.e., you must install the certificate within 14 days from the time you requested it. cobas link Configurator A component of the cobas link software. Used by Roche Service representatives and the administrator to set up and to administrate a cobas link, e.g., to download patches manually. 13 cobas link Essential Information Upload Maintenance item that must be included in the daily wake-up pipe. Performs a backup of the instrument master data to the cobas link data station. cobas link stores the last five backups, current reagent Reagent packs that are currently in use, also called active reagent packs. e-barcode Contains machine-readable information for cobas instruments: test-specific instrument settings, lot-specific instrument settings, and values of calibrators and controls. e-library package Information unit downloaded from the remote service platform to a cobas link. Contains e-library package inserts and e-barcodes. e-library package inserts are for the operator. e-barcodes contain machine-readable parameters that are passed on to the instruments.

e-library package insert The following documents come as e-library package inserts: instructions for use, value sheets for calibrators and controls, important notes, and announcements from the local technical support. E.Stop Analyzer status displayed in the status line. The analyzer performed an emergency stop due to hardware failure or because any of the safety devices requested an emergency stop. The analyzer requires either a complete power-off or a reset to resume operation. global button Feature of the user interface. Global buttons are displayed with all menus and provide quick access to key functions of the analyzer. HHT driver Hitachi driver. A component of the cobas link software. Manages the information flow between the cobas link and the cobas instruments, e.g., the release of e-barcodes to the instruments. HRD viewer Human-Readable Data viewer. A component of the cobas link software. Provides the user interface of the cobas e-library application. LIS Laboratory Information System. Passes orders to the instruments and receives test results. maintenance mode Analyzer mode activated for doing maintenance work. The maintenance switch is on. You may open the top cover as all components on top are switched off.

The status line is yellow and displays the Standby status. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 370 Glossary online help - System Overview menu online help Provides context-sensitive access via F1 to the Operator's Manual and the Software Description . Pause/Scan button Global button used to pause sampling, e.g., to insert a STAT sample, or to scan the sample containers on the sample disk. Photometer Check report Report generated by the Photometer check maintenance item. The photometer report displays the current and previous absorbance values for all wavelengths. The report may indicate a problem with the reaction cells, the photometric window, or the photometer lamp. Make sure that absorbance values do not exceed 14 000 and the current values do not differ significantly from previous results. RCL Roche Connectivity Layer. A component of the cobas link software. Manages the information flow between the cobas link and the remote service platform. Remote service platform Provides information services for operators and instruments via a secure Internet connection. Services include the following: • Downloads of the most recent information on applications, calibrators, and controls • Support and troubleshooting via screen sharing • Advanced performance monitoring by evaluating uploads of statistical data from the instruments S.Stop button Global button used to stop sampling. The analyzer completes measuring already pipetted samples. Then, the analyzer switches to Standby status. Shutdown button Global button used to log off, start Sleep mode, or shut down the analyzer.

standby reagent Reagent packs that are already on board but not yet in use. status line Feature of the user interface. Displays the analyzer status, login name, date/time, and time remaining to finish the current maintenance task. System Overview menu Root menu of the user interface. Used to prepare the analyzer for daily operation. Displays the button about the modules of the analyzer and provides access to the five main menus. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Index 371 Index Symbols B (10) ISE Wash, 224, 261 (12) Flow Path Wash, 289 (2) Check, 284 (3) Cell Blank Measurement, 276, 319 (3) Mechanisms Check, 345, 350 (4) Incubation Water Exchange, 282, 302 (5) Air Purge, 335, 344 (6) Cuvette Mixing, 284 (6) Wash Reaction Parts, 275 (9) Incubator Bath Cleaning, 302 Numerics 2 Point, calibration method, 122 A Abbreviations, list of, 10 Accessory - printer, 21 ACN (Application Code Number), 203 Air purge - check, 66 - performing, 335, 344 Alarm - check system alarms, 75 - defining reagent level alarms, 116 - volume knob, 22 Alarm button - flashing, 75 Analyzer - checking, 66, 72 - cleaning surface, 268 - dimensions, 359 - preparing for operation, 76 - quick start guide, 58 - shutting down, 227 - specifications, 357 - starting from power-off, 70 - statuses, 39 - switching off, 221, 223 - switching on, 70 - waking up, 63 Analyzer, areas of - ISE area, 19 - reaction area, 19 - reagent area, 19 - sampling area, 19 Application - downloading new application parameters, 86 Application Code Number (ACN), 203 Approvals, 4 Backing up data, 78 Backup - on USB flash drive, 190 - patient and QC results, 191 Backup Data button, 191 Backup Data dialog box, 78 Barcode - barcoded and non-barcoded samples, 212 - labels, 366 - placing labels on sample tubes, 213 - reagent packs, 102 - specifications, list of, 366 - unreadable, 213 Barcode read errors, managing, 213 Blank, calibration method, 122 Buttons - alarm button, 75 - global buttons, 32 - help button, 42 - shapes, 38 C Calibration - checking results, 132 - introduction, 121 - performing manual calibration, 126, 129 - performing recommended calibration, 124, 127 - requesting, 82, 127 Calibration and QC - performing, 81 - requesting, 82 Calibration and QC Load List, 83 Calibration And QC Select button, 82, 124, 127, 144 Calibration menu - Calibrator submenu, 133 - Status submenu, 129, 132 Calibration methods, 122 Calibrator - assigning non-barcoded calibrators, 133 - overview, 96 Calibrators and controls - measuring, 84 - printing load list, 83 Cell blank measurement - performing, 276, 319 Cell cleaners - Cell cleaner I, 22 Cell covers - cleaning, 270 Cell detergent prime, performing, 115 Cell detergents cobas link, e.g., to download patches manually, 13 cobas link Essential Information Upload Maintenance item that must be included in the daily wake-up pipe. Performs a backup of the instrument master data to the cobas link data station. cobas link stores the last five backups, current reagent Reagent packs that are currently in use, also called active reagent packs. e-barcode Contains machine-readable information for cobas instruments: test-specific instrument settings, lot-specific instrument settings, and values of calibrators and controls. e-library package Information unit downloaded from the remote service platform to a cobas link. Contains e-library package inserts and e-barcodes. e-library package inserts are for the operator. e-barcodes contain machine-readable parameters that are passed on to the instruments.

e-library package insert The following documents come as e-library package inserts: instructions for use, value sheets for calibrators and controls, important notes, and announcements from the local technical support. E.Stop Analyzer status displayed in the status line. The analyzer performed an emergency stop due to hardware failure or because any of the safety devices requested an emergency stop. The analyzer requires either a complete power-off or a reset to resume operation. global button Feature of the user interface. Global buttons are displayed with all menus and provide quick access to key functions of the analyzer. HHT driver Hitachi driver. A component of the cobas link software. Manages the information flow between the cobas link and the cobas instruments, e.g., the release of e-barcodes to the instruments. HRD viewer Human-Readable Data viewer. A component of the cobas link software. Provides the user interface of the cobas e-library application. LIS Laboratory Information System. Passes orders to the instruments and receives test results. maintenance mode Analyzer mode activated for doing maintenance work. The maintenance switch is on. You may open the top cover as all components on top are switched off.

The status line is yellow and displays the Standby status. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 370 Glossary online help - System Overview menu online help Provides context-sensitive access via F1 to the Operator's Manual and the Software Description . Pause/Scan button Global button used to pause sampling, e.g., to insert a STAT sample, or to scan the sample containers on the sample disk. Photometer Check report Report generated by the Photometer check maintenance item. The photometer report displays the current and previous absorbance values for all wavelengths. The report may indicate a problem with the reaction cells, the photometric window, or the photometer lamp. Make sure that absorbance values do not exceed 14 000 and the current values do not differ significantly from previous results. RCL Roche Connectivity Layer. A component of the cobas link software. Manages the information flow between the cobas link and the remote service platform. Remote service platform Provides information services for operators and instruments via a secure Internet connection. Services include the following: • Downloads of the most recent information on applications, calibrators, and controls • Support and troubleshooting via screen sharing • Advanced performance monitoring by evaluating uploads of statistical data from the instruments S.Stop button Global button used to stop sampling. The analyzer completes measuring already pipetted samples. Then, the analyzer switches to Standby status. Shutdown button Global button used to log off, start Sleep mode, or shut down the analyzer.

standby reagent Reagent packs that are already on board but not yet in use. status line Feature of the user interface. Displays the analyzer status, login name, date/time, and time remaining to finish the current maintenance task. System Overview menu Root menu of the user interface. Used to prepare the analyzer for daily operation. Displays the button about the modules of the analyzer and provides access to the five main menus. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Index 371 Index Symbols B (10) ISE Wash, 224, 261 (12) Flow Path Wash, 289 (2) Check, 284 (3) Cell Blank Measurement, 276, 319 (3) Mechanisms Check, 345, 350 (4) Incubation Water Exchange, 282, 302 (5) Air Purge, 335, 344 (6) Cuvette Mixing, 284 (6) Wash Reaction Parts, 275 (9) Incubator Bath Cleaning, 302 Numerics 2 Point, calibration method, 122 A Abbreviations, list of, 10 Accessory - printer, 21 ACN (Application Code Number), 203 Air purge - check, 66 - performing, 335, 344 Alarm - check system alarms, 75 - defining reagent level alarms, 116 - volume knob, 22 Alarm button - flashing, 75 Analyzer - checking, 66, 72 - cleaning surface, 268 - dimensions, 359 - preparing for operation, 76 - quick start guide, 58 - shutting down, 227 - specifications, 357 - starting from power-off, 70 - statuses, 39 - switching off, 221, 223 - switching on, 70 - waking up, 63 Analyzer, areas of - ISE area, 19 - reaction area, 19 - reagent area, 19 - sampling area, 19 Application - downloading new application parameters, 86 Application Code Number (ACN), 203 Approvals, 4 Backing up data, 78 Backup - on USB flash drive, 190 - patient and QC results, 191 Backup Data button, 191 Backup Data dialog box, 78 Barcode - barcoded and non-barcoded samples, 212 - labels, 366 - placing labels on sample tubes, 213 - reagent packs, 102 - specifications, list of, 366 - unreadable, 213 Barcode read errors, managing, 213 Blank, calibration method, 122 Buttons - alarm button, 75 - global buttons, 32 - help button, 42 - shapes, 38 C Calibration - checking results, 132 - introduction, 121 - performing manual calibration, 126, 129 - performing recommended calibration, 124, 127 - requesting, 82, 127 Calibration and QC - performing, 81 - requesting, 82 Calibration and QC Load List, 83 Calibration And QC Select button, 82, 124, 127, 144 Calibration menu - Calibrator submenu, 133 - Status submenu, 129, 132 Calibration methods, 122 Calibrator - assigning non-barcoded calibrators, 133 - overview, 96 Calibrators and controls - measuring, 84 - printing load list, 83 Cell blank measurement - performing, 276, 319 Cell cleaners - Cell cleaner I, 22 Cell covers - cleaning, 270 Cell detergent prime, performing, 115 Cell detergents cobas link, e.g., to download patches manually, 13 cobas link Essential Information Upload Maintenance item that must be included in the daily wake-up pipe. Performs a backup of the instrument master data to the cobas link data station. cobas link stores the last five backups, current reagent Reagent packs that are currently in use, also called active reagent packs. e-barcode Contains machine-readable information for cobas instruments: test-specific instrument settings, lot-specific instrument settings, and values of calibrators and controls. e-library package Information unit downloaded from the remote service platform to a cobas link. Contains e-library package inserts and e-barcodes. e-library package inserts are for the operator. e-barcodes contain machine-readable parameters that are passed on to the instruments.

e-library package insert The following documents come as e-library package inserts: instructions for use, value sheets for calibrators and controls, important notes, and announcements from the local technical support. E.Stop Analyzer status displayed in the status line. The analyzer performed an emergency stop due to hardware failure or because any of the safety devices requested an emergency stop. The analyzer requires either a complete power-off or a reset to resume operation. global button Feature of the user interface. Global buttons are displayed with all menus and provide quick access to key functions of the analyzer. HHT driver Hitachi driver. A component of the cobas link software. Manages the information flow between the cobas link and the cobas instruments, e.g., the release of e-barcodes to the instruments. HRD viewer Human-Readable Data viewer. A component of the cobas link software. Provides the user interface of the cobas e-library application. LIS Laboratory Information System. Passes orders to the instruments and receives test results. maintenance mode Analyzer mode activated for doing maintenance work. The maintenance switch is on. You may open the top cover as all components on top are switched off.

The status line is yellow and displays the Standby status. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 370 Glossary online help - System Overview menu online help Provides context-sensitive access via F1 to the Operator's Manual and the Software Description . Pause/Scan button Global button used to pause sampling, e.g., to insert a STAT sample, or to scan the sample containers on the sample disk. Photometer Check report Report generated by the Photometer check maintenance item. The photometer report displays the current and previous absorbance values for all wavelengths. The report may indicate a problem with the reaction cells, the photometric window, or the photometer lamp. Make sure that absorbance values do not exceed 14 000 and the current values do not differ significantly from previous results. RCL Roche Connectivity Layer. A component of the cobas link software. Manages the information flow between the cobas link and the remote service platform. Remote service platform Provides information services for operators and instruments via a secure Internet connection. Services include the following: • Downloads of the most recent information on applications, calibrators, and controls • Support and troubleshooting via screen sharing • Advanced performance monitoring by evaluating uploads of statistical data from the instruments S.Stop button Global button used to stop sampling. The analyzer completes measuring already pipetted samples. Then, the analyzer switches to Standby status. Shutdown button Global button used to log off, start Sleep mode, or shut down the analyzer.

standby reagent Reagent packs that are already on board but not yet in use. status line Feature of the user interface. Displays the analyzer status, login name, date/time, and time remaining to finish the current maintenance task. System Overview menu Root menu of the user interface. Used to prepare the analyzer for daily operation. Displays the button about the modules of the analyzer and provides access to the five main menus. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Index 371 Index Symbols B (10) ISE Wash, 224, 261 (12) Flow Path Wash, 289 (2) Check, 284 (3) Cell Blank Measurement, 276, 319 (3) Mechanisms Check, 345, 350 (4) Incubation Water Exchange, 282, 302 (5) Air Purge, 335, 344 (6) Cuvette Mixing, 284 (6) Wash Reaction Parts, 275 (9) Incubator Bath Cleaning, 302 Numerics 2 Point, calibration method, 122 A Abbreviations, list of, 10 Accessory - printer, 21 ACN (Application Code Number), 203 Air purge - check, 66 - performing, 335, 344 Alarm - check system alarms, 75 - defining reagent level alarms, 116 - volume knob, 22 Alarm button - flashing, 75 Analyzer - checking, 66, 72 - cleaning surface, 268 - dimensions, 359 - preparing for operation, 76 - quick start guide, 58 - shutting down, 227 - specifications, 357 - starting from power-off, 70 - statuses, 39 - switching off, 221, 223 - switching on, 70 - waking up, 63 Analyzer, areas of - ISE area, 19 - reaction area, 19 - reagent area, 19 - sampling area, 19 Application - downloading new application parameters, 86 Application Code Number (ACN), 203 Approvals, 4 Backing up data, 78 Backup - on USB flash drive, 190 - patient and QC results, 191 Backup Data button, 191 Backup Data dialog box, 78 Barcode - barcoded and non-barcoded samples, 212 - labels, 366 - placing labels on sample tubes, 213 - reagent packs, 102 - specifications, list of, 366 - unreadable, 213 Barcode read errors, managing, 213 Blank, calibration method, 122 Buttons - alarm button, 75 - global buttons, 32 - help button, 42 - shapes, 38 C Calibration - checking results, 132 - introduction, 121 - performing manual calibration, 126, 129 - performing recommended calibration, 124, 127 - requesting, 82, 127 Calibration and QC - performing, 81 - requesting, 82 Calibration and QC Load List, 83 Calibration And QC Select button, 82, 124, 127, 144 Calibration menu - Calibrator submenu, 133 - Status submenu, 129, 132 Calibration methods, 122 Calibrator - assigning non-barcoded calibrators, 133 - overview, 96 Calibrators and controls - measuring, 84 - printing load list, 83 Cell blank measurement - performing, 276, 319 Cell cleaners - Cell cleaner I, 22 Cell covers - cleaning, 270 Cell detergent prime, performing, 115 Cell detergents cobas link, e.g., to download patches manually, 13 cobas link Essential Information Upload Maintenance item that must be included in the daily wake-up pipe. Performs a backup of the instrument master data to the cobas link data station. cobas link stores the last five backups, current reagent Reagent packs that are currently in use, also called active reagent packs. e-barcode Contains machine-readable information for cobas instruments: test-specific instrument settings, lot-specific instrument settings, and values of calibrators and controls. e-library package Information unit downloaded from the remote service platform to a cobas link. Contains e-library package inserts and e-barcodes. e-library package inserts are for the operator. e-barcodes contain machine-readable parameters that are passed on to the instruments.

e-library package insert The following documents come as e-library package inserts: instructions for use, value sheets for calibrators and controls, important notes, and announcements from the local technical support. E.Stop Analyzer status displayed in the status line. The analyzer performed an emergency stop due to hardware failure or because any of the safety devices requested an emergency stop. The analyzer requires either a complete power-off or a reset to resume operation. global button Feature of the user interface. Global buttons are displayed with all menus and provide quick access to key functions of the analyzer. HHT driver Hitachi driver. A component of the cobas link software. Manages the information flow between the cobas link and the cobas instruments, e.g., the release of e-barcodes to the instruments. HRD viewer Human-Readable Data viewer. A component of the cobas link software. Provides the user interface of the cobas e-library application. LIS Laboratory Information System. Passes orders to the instruments and receives test results. maintenance mode Analyzer mode activated for doing maintenance work. The maintenance switch is on. You may open the top cover as all components on top are switched off.

The status line is yellow and displays the Standby status. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 370 Glossary online help - System Overview menu online help Provides context-sensitive access via F1 to the Operator's Manual and the Software Description . Pause/Scan button Global button used to pause sampling, e.g., to insert a STAT sample, or to scan the sample containers on the sample disk. Photometer Check report Report generated by the Photometer check maintenance item. The photometer report displays the current and previous absorbance values for all wavelengths. The report may indicate a problem with the reaction cells, the photometric window, or the photometer lamp. Make sure that absorbance values do not exceed 14 000 and the current values do not differ significantly from previous results. RCL Roche Connectivity Layer. A component of the cobas link software. Manages the information flow between the cobas link and the remote service platform. Remote service platform Provides information services for operators and instruments via a secure Internet connection. Services include the following: • Downloads of the most recent information on applications, calibrators, and controls • Support and troubleshooting via screen sharing • Advanced performance monitoring by evaluating uploads of statistical data from the instruments S.Stop button Global button used to stop sampling. The analyzer completes measuring already pipetted samples. Then, the analyzer switches to Standby status. Shutdown button Global button used to log off, start Sleep mode, or shut down the analyzer.

standby reagent Reagent packs that are already on board but not yet in use. status line Feature of the user interface. Displays the analyzer status, login name, date/time, and time remaining to finish the current maintenance task. System Overview menu Root menu of the user interface. Used to prepare the analyzer for daily operation. Displays the button about the modules of the analyzer and provides access to the five main menus. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Index 371 Index Symbols B (10) ISE Wash, 224, 261 (12) Flow Path Wash, 289 (2) Check, 284 (3) Cell Blank Measurement, 276, 319 (3) Mechanisms Check, 345, 350 (4) Incubation Water Exchange, 282, 302 (5) Air Purge, 335, 344 (6) Cuvette Mixing, 284 (6) Wash Reaction Parts, 275 (9) Incubator Bath Cleaning, 302 Numerics 2 Point, calibration method, 122 A Abbreviations, list of, 10 Accessory - printer, 21 ACN (Application Code Number), 203 Air purge - check, 66 - performing, 335, 344 Alarm - check system alarms, 75 - defining reagent level alarms, 116 - volume knob, 22 Alarm button - flashing, 75 Analyzer - checking, 66, 72 - cleaning surface, 268 - dimensions, 359 - preparing for operation, 76 - quick start guide, 58 - shutting down, 227 - specifications, 357 - starting from power-off, 70 - statuses, 39 - switching off, 221, 223 - switching on, 70 - waking up, 63 Analyzer, areas of - ISE area, 19 - reaction area, 19 - reagent area, 19 - sampling area, 19 Application - downloading new application parameters, 86 Application Code Number (ACN), 203 Approvals, 4 Backing up data, 78 Backup - on USB flash drive, 190 - patient and QC results, 191 Backup Data button, 191 Backup Data dialog box, 78 Barcode - barcoded and non-barcoded samples, 212 - labels, 366 - placing labels on sample tubes, 213 - reagent packs, 102 - specifications, list of, 366 - unreadable, 213 Barcode read errors, managing, 213 Blank, calibration method, 122 Buttons - alarm button, 75 - global buttons, 32 - help button, 42 - shapes, 38 C Calibration - checking results, 132 - introduction, 121 - performing manual calibration, 126, 129 - performing recommended calibration, 124, 127 - requesting, 82, 127 Calibration and QC - performing, 81 - requesting, 82 Calibration and QC Load List, 83 Calibration And QC Select button, 82, 124, 127, 144 Calibration menu - Calibrator submenu, 133 - Status submenu, 129, 132 Calibration methods, 122 Calibrator - assigning non-barcoded calibrators, 133 - overview, 96 Calibrators and controls - measuring, 84 - printing load list, 83 Cell blank measurement - performing, 276, 319 Cell cleaners - Cell cleaner I, 22 Cell covers - cleaning, 270 Cell detergent prime, performing, 115 Cell detergents cobas link, e.g., to download patches manually, 13 cobas link Essential Information Upload Maintenance item that must be included in the daily wake-up pipe. Performs a backup of the instrument master data to the cobas link data station. cobas link stores the last five backups, current reagent Reagent packs that are currently in use, also called active reagent packs. e-barcode Contains machine-readable information for cobas instruments: test-specific instrument settings, lot-specific instrument settings, and values of calibrators and controls. e-library package Information unit downloaded from the remote service platform to a cobas link. Contains e-library package inserts and e-barcodes. e-library package inserts are for the operator. e-barcodes contain machine-readable parameters that are passed on to the instruments.

e-library package insert The following documents come as e-library package inserts: instructions for use, value sheets for calibrators and controls, important notes, and announcements from the local technical support. E.Stop Analyzer status displayed in the status line. The analyzer performed an emergency stop due to hardware failure or because any of the safety devices requested an emergency stop. The analyzer requires either a complete power-off or a reset to resume operation. global button Feature of the user interface. Global buttons are displayed with all menus and provide quick access to key functions of the analyzer. HHT driver Hitachi driver. A component of the cobas link software. Manages the information flow between the cobas link and the cobas instruments, e.g., the release of e-barcodes to the instruments. HRD viewer Human-Readable Data viewer. A component of the cobas link software. Provides the user interface of the cobas e-library application. LIS Laboratory Information System. Passes orders to the instruments and receives test results. maintenance mode Analyzer mode activated for doing maintenance work. The maintenance switch is on. You may open the top cover as all components on top are switched off.

The status line is yellow and displays the Standby status. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 370 Glossary online help - System Overview menu online help Provides context-sensitive access via F1 to the Operator's Manual and the Software Description . Pause/Scan button Global button used to pause sampling, e.g., to insert a STAT sample, or to scan the sample containers on the sample disk. Photometer Check report Report generated by the Photometer check maintenance item. The photometer report displays the current and previous absorbance values for all wavelengths. The report may indicate a problem with the reaction cells, the photometric window, or the photometer lamp. Make sure that absorbance values do not exceed 14 000 and the current values do not differ significantly from previous results. RCL Roche Connectivity Layer. A component of the cobas link software. Manages the information flow between the cobas link and the remote service platform. Remote service platform Provides information services for operators and instruments via a secure Internet connection. Services include the following: • Downloads of the most recent information on applications, calibrators, and controls • Support and troubleshooting via screen sharing • Advanced performance monitoring by evaluating uploads of statistical data from the instruments S.Stop button Global button used to stop sampling. The analyzer completes measuring already pipetted samples. Then, the analyzer switches to Standby status. Shutdown button Global button used to log off, start Sleep mode, or shut down the analyzer.

standby reagent Reagent packs that are already on board but not yet in use. status line Feature of the user interface. Displays the analyzer status, login name, date/time, and time remaining to finish the current maintenance task. System Overview menu Root menu of the user interface. Used to prepare the analyzer for daily operation. Displays the button about the modules of the analyzer and provides access to the five main menus. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Index 371 Index Symbols B (10) ISE Wash, 224, 261 (12) Flow Path Wash, 289 (2) Check, 284 (3) Cell Blank Measurement, 276, 319 (3) Mechanisms Check, 345, 350 (4) Incubation Water Exchange, 282, 302 (5) Air Purge, 335, 344 (6) Cuvette Mixing, 284 (6) Wash Reaction Parts, 275 (9) Incubator Bath Cleaning, 302 Numerics 2 Point, calibration method, 122 A Abbreviations, list of, 10 Accessory - printer, 21 ACN (Application Code Number), 203 Air purge - check, 66 - performing, 335, 344 Alarm - check system alarms, 75 - defining reagent level alarms, 116 - volume knob, 22 Alarm button - flashing, 75 Analyzer - checking, 66, 72 - cleaning surface, 268 - dimensions, 359 - preparing for operation, 76 - quick start guide, 58 - shutting down, 227 - specifications, 357 - starting from power-off, 70 - statuses, 39 - switching off, 221, 223 - switching on, 70 - waking up, 63 Analyzer, areas of - ISE area, 19 - reaction area, 19 - reagent area, 19 - sampling area, 19 Application - downloading new application parameters, 86 Application Code Number (ACN), 203 Approvals, 4 Backing up data, 78 Backup - on USB flash drive, 190 - patient and QC results, 191 Backup Data button, 191 Backup Data dialog box, 78 Barcode - barcoded and non-barcoded samples, 212 - labels, 366 - placing labels on sample tubes, 213 - reagent packs, 102 - specifications, list of, 366 - unreadable, 213 Barcode read errors, managing, 213 Blank, calibration method, 122 Buttons - alarm button, 75 - global buttons, 32 - help button, 42 - shapes, 38 C Calibration - checking results, 132 - introduction, 121 - performing manual calibration, 126, 129 - performing recommended calibration, 124, 127 - requesting, 82, 127 Calibration and QC - performing, 81 - requesting, 82 Calibration and QC Load List, 83 Calibration And QC Select button, 82, 124, 127, 144 Calibration menu - Calibrator submenu, 133 - Status submenu, 129, 132 Calibration methods, 122 Calibrator - assigning non-barcoded calibrators, 133 - overview, 96 Calibrators and controls - measuring, 84 - printing load list, 83 Cell blank measurement - performing, 276, 319 Cell cleaners - Cell cleaner I, 22 Cell covers - cleaning, 270 Cell detergent prime, performing, 115 Cell detergents cobas link, e.g., to download patches manually, 13 cobas link Essential Information Upload Maintenance item that must be included in the daily wake-up pipe. Performs a backup of the instrument master data to the cobas link data station. cobas link stores the last five backups, current reagent Reagent packs that are currently in use, also called active reagent packs. e-barcode Contains machine-readable information for cobas instruments: test-specific instrument settings, lot-specific instrument settings, and values of calibrators and controls. e-library package Information unit downloaded from the remote service platform to a cobas link. Contains e-library package inserts and e-barcodes. e-library package inserts are for the operator. e-barcodes contain machine-readable parameters that are passed on to the instruments.

e-library package insert The following documents come as e-library package inserts: instructions for use, value sheets for calibrators and controls, important notes, and announcements from the local technical support. E.Stop Analyzer status displayed in the status line. The analyzer performed an emergency stop due to hardware failure or because any of the safety devices requested an emergency stop. The analyzer requires either a complete power-off or a reset to resume operation. global button Feature of the user interface. Global buttons are displayed with all menus and provide quick access to key functions of the analyzer. HHT driver Hitachi driver. A component of the cobas link software. Manages the information flow between the cobas link and the cobas instruments, e.g., the release of e-barcodes to the instruments. HRD viewer Human-Readable Data viewer. A component of the cobas link software. Provides the user interface of the cobas e-library application. LIS Laboratory Information System. Passes orders to the instruments and receives test results. maintenance mode Analyzer mode activated for doing maintenance work. The maintenance switch is on. You may open the top cover as all components on top are switched off.

The status line is yellow and displays the Standby status. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 370 Glossary online help - System Overview menu online help Provides context-sensitive access via F1 to the Operator's Manual and the Software Description . Pause/Scan button Global button used to pause sampling, e.g., to insert a STAT sample, or to scan the sample containers on the sample disk. Photometer Check report Report generated by the Photometer check maintenance item. The photometer report displays the current and previous absorbance values for all wavelengths. The report may indicate a problem with the reaction cells, the photometric window, or the photometer lamp. Make sure that absorbance values do not exceed 14 000 and the current values do not differ significantly from previous results. RCL Roche Connectivity Layer. A component of the cobas link software. Manages the information flow between the cobas link and the remote service platform. Remote service platform Provides information services for operators and instruments via a secure Internet connection. Services include the following: • Downloads of the most recent information on applications, calibrators, and controls • Support and troubleshooting via screen sharing • Advanced performance monitoring by evaluating uploads of statistical data from the instruments S.Stop button Global button used to stop sampling. The analyzer completes measuring already pipetted samples. Then, the analyzer switches to Standby status. Shutdown button Global button used to log off, start Sleep mode, or shut down the analyzer.

standby reagent Reagent packs that are already on board but not yet in use. status line Feature of the user interface. Displays the analyzer status, login name, date/time, and time remaining to finish the current maintenance task. System Overview menu Root menu of the user interface. Used to prepare the analyzer for daily operation. Displays the button about the modules of the analyzer and provides access to the five main menus. Roche Diagnostics cobas c 311 analyzer · Software version 01-10 · Operator's Manual · 3.1 Index 371 Index Symbols B (10) ISE Wash, 224, 261 (12) Flow Path Wash, 289 (2) Check, 284 (3) Cell Blank Measurement, 276, 319 (3) Mechanisms Check, 345, 350 (4) Incubation Water Exchange, 282, 302 (5) Air Purge, 335, 344 (6) Cuvette Mixing, 284 (6) Wash Reaction Parts, 275 (9) Incubator Bath Cleaning, 302 Numerics 2 Point, calibration method, 122 A Abbreviations, list of, 10 Accessory - printer, 21 ACN (Application Code Number), 203 Air purge - check, 66 - performing, 335, 344 Alarm - check system alarms, 75 - defining reagent level alarms, 116 - volume knob, 22 Alarm button - flashing, 75 Analyzer - checking, 66, 72 - cleaning surface, 268 - dimensions, 359 - preparing for operation, 76 - quick start guide, 58 - shutting down, 227 - specifications, 357 - starting from power-off, 70 - statuses, 39 - switching off, 221, 223 - switching on, 70 - waking up, 63 Analyzer, areas of - ISE area, 19 - reaction area, 19 - reagent area, 19 - sampling area, 19 Application - downloading new application parameters, 86 Application Code Number (ACN), 203 Approvals, 4 Backing up data, 78 Backup - on USB flash drive, 190 - patient and QC results, 191 Backup Data button, 191 Backup Data dialog box